Annex 1:

Standard Operation Procedures (SOP)

For Environmental Monitoring

1.2 Chemical and Biological Water Quality (1)

Standard Operation Procedure (SOP) For the Determination of Oil and Grease

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Prepared by: _____ Date: _____
Chemist

Reviewed by: _____ Date: _____
Laboratory Chief or Quality Assurance Coordinator

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Collect a representative grab sample in a wide-mouth glass bottle that has been washed with soap, rinsed with water, and finally rinsed with solvent to remove any residues that might interfere with the analysis. Do not overfill the sample container and do not subdivide the sample in the laboratory.

Typically, collect wastewater samples of approximately 1 L. If sample concentration is expected to be greater than 1,000 mg extractable material/L, collect proportionately smaller volumes. If analysis is to be delayed for more than 2 hr, acidify to pH 2 or lower with either 1:1 HCl or 1:1 H₂SO₄ and refrigerate.

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

i.e. the quantity of the used sample equals two folds of the used solvent (where C is the oil concentration in the sample) which is transformed to the solvent phase and the concentration becomes 2C (which is corresponding to the standard solution we have prepared), and required concentration (span value) during span calibration is half; i.e 100 mg/L

→ when the standard solution concentration is 200 mg/L, the corresponding span should be 100

3- Procedures

Repeat stages of zero calibration as follows:
 take 10ml of solvent in a 10ml syringe after washing using the same solvent 3 times, and inject it to the measurement chamber
 - take 20ml of distilled or de-ionized water in a 20ml syringe after washing it 3 times with the same water, and inject it to the measurement chamber
 - add few drops of 1:1 hydrochloric acid (HCl) to accelerate the separation process
 - press [EXTRACT] on the panel, and wait until timer annulated (usually 40 Sec.), and then wait 10-15 Sec.
 - open the top valve to the right and wait 10-15 Sec. then press [MEASURE] button until the reading is stabilized, and when the green color is lit press [HOLD] then continue as following:

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

Taking Measurements without using the extraction units

Notice: It is important to remember that solvent extraction is conducted outside of the device, therefore we don't use the extraction function of the device, instead, we conduct direct reading. This method is applied when the sample contains big amount of oil and floating grease which could be noticed by bare eye, and accompanied with a lot of turbidity.


Necessary tools and chemicals

- Separation funnel (300-500 ml)
- Cylinder (100 ml, 200 ml)
- 1:1 Hydrochloric acid (HCl)
- Solvent (S-316)
- PH meter
- Other materials (tissue paper, non-aqueous sodium-soleplate or equivalent as required)

Repeat this operation more than 3 times, as obviously, these steps are an iteration of the basic measurement procedures.


$$\text{Result} = \text{Average of Measurements} \times \text{Dilution factor}$$

Note You cannot obtain a correct calibration if the parts that come into contact with the liquid are contaminated with liquid remaining from a previous operation. Repeat the liquid transfer → drainage (purge) cycle at least 3 times, before executing the liquid transfer → measurement → drainage (measurement).

- Beakers, Syringes, Cylinders, and all glassware used for solvent, should be rinsed with solvent.
- Syringe used for water, should be washed with de-ionized water

● FLOW DIAGRAM

No.	Parts name
1	Sample inlet
2	Stirrer
3	Extracting chamber
4	Filter
5	Ball valve
6	Drain
7	Air Vent
8	Analyzer unit

Short Title: Phosphorus, Reactive (Orthophosphate)

Revision No.: 3

Date: Aug. 2007

Page 1 of 8

**Standard Operation Procedure (SOP)
For the Determination of
Phosphorus, Reactive (Orthophosphate)**

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Prepared by: _____ Date:
Chemist

Reviewed by: _____ Date:
Laboratory Chief or Quality Assurance Coordinator

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Date: Aug. 2007

Page 3 of 8

containing phosphate for cleaning glassware used in phosphate analysis.

Analyze samples immediately for best results. If prompt analysis is not possible, preserve samples by filtering immediately and storing at 4°C for up to 48 hours. The sample should be at room temperature before analysis.

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1 - Reagent blank

For more accurate results, determine a reagent blank value for each new lot of reagent. Follow the procedure using deionized water instead of the sample. Subtract the reagent blank value from the final results or perform a reagent blank adjust.

2- Standard Solution Method

1. Prepare a 2.00 mg/L phosphate standard by pipetting 4.00 mL of 50 mg/L Phosphate Standard Solution into a 100 mL volumetric flask. Dilute to volume with demineralized water and mix. Use this solution in place of the sample, and perform the test as described above. (Alternately, use one of the mixed parameter standards listed in. These contain 2.0 mg/L phosphate.)

(Preparation of Standard phosphate solution: Dissolve in distilled water 71.64 mg anhydrous KH_2PO_4 and dilute to 1,000 mL; 1.00 mL = 50.0 $\mu g PO_4^{3-}$, 50 mg PO_4^{3-}/L)

2. To adjust the calibration curve using the reading obtained with the standard solution, press **OPTIONS>MORE** on the current program menu. Press **STANDARD ADJUST**.

3. Press **ON**. Press **ADJUST** to accept the displayed concentration. If an alternate concentration is used, press the number in the box to enter the actual concentration, then press **OK**. Press **ADJUST**.

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1. After reading test results, leave the sample cell (unspiked sample) in the instrument. Verify the chemical form.
2. Press **OPTIONS>MORE**. Press **STANDARD ADDITIONS**. A summary of the standard additions procedure will appear.
3. Press **OK** to accept the default values for standard concentration, sample volume, and spike volumes. Press **EDIT** to change these values. After values are accepted, the unspiked sample reading will appear in the top row.
4. Open a Phosphate 10-mL Ampule Standard, 50-mg/L po₄³⁻
5. Prepare a 0.1-mL sample spike by adding 0.1 mL of standard to the unspiked sample. Press the timer icon. After the timer beeps, read the result.
6. Prepare a 0.2-mL sample spike by adding 0.1 mL of standard to the 0.1-mL sample spike. Press the timer icon. After the timer beeps, read the result.

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7. Prepare a 0.3-mL sample spike by adding 0.1 mL of standard to the 0.2-mL sample spike.

Press the timer icon. After the timer beeps, read the result. Each addition should reflect approximately 100% recovery.

Note: For AccuVac® ampules, fill three Mixing Cylinders† with 50 mL of sample and spike with 0.2 mL, 0.4 mL, and 0.6 mL of standard. Transfer 40 mL from each of the three mixing cylinders to three 50-mL beakers*. Analyze each standard addition sample as described in the procedure above. Accept each standard additions reading by pressing **READ**. Each addition should reflect approximately 100% recovery.

8. After completing the sequence, press **GRAPH** to view the best-fit line through the standard additions data points, accounting for matrix interferences. Press **IDEAL LINE** to view Relationships between the sample spikes and the “Ideal Line” of 100% recovery.

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


**Standard Operation Procedure (SOP)
For the Determination of Chromium, Total**

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1. Scope and Application: water, wastewater

2. Summary of Method: Alkaline Hypobromite Oxidation method

3. Measurement range: 0.01 to 0.70 mg/L

4. Necessary Equipment and Supplies:

- 1) UV/VIS Spectrophotometer (HACH DR 5000)
- 2) Sample Cells, 1-inch square, 10mL, matched pair: 2/Test
- 3) Sample Cell, 10-20-25 mL, with cap: 1
- 4) Hot plate, Water bath and Rack: 1

5. Required Reagents

- 1) Acid Reagent Powder Pillows (Cat No.: 2126-99): 1/Test
- 2) Chroma Ver[®]3 Chromium Reagent Powder Pillows (Cat. No. 12066-99): 1/Test
- 3) Chromium 1 Reagent Powder Pillows (Cat. No.: 2043-99): 1/Test
- 4) Chromium 2 Reagent Powder Pillows (Cat. No.: 2044-99): 1/Test

6. Principle of Determination Trivalent chromium in the sample is oxidized to the Hexavalent form by hypobromite ion under alkaline conditions. The sample is acidified. The total chromium content is determined by the 1,5-Diphenylcarbohydrazide method. Determine trivalent chromium by subtracting the results of a separate Hexavalent chromium test from the results of the total chromium test. Test results are measured at 540 nm.

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7. Sample Collection, Storage, and Preservation

Collect sample in acid-washed glass or plastic containers. To preserve samples, adjust the pH to 2 or less with nitric acid. This requires approximately 2 mL per liter of the acid. Store preserved samples at room temperature up to six months.

Adjust the pH to about 4 with 5.0 N Sodium Hydroxide before analysis. Correct the test result for volume additions.

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1 - Reagent blank

For more accurate results, determine a reagent blank value for each new lot of reagent. Follow the procedure using deionized water in place of the sample. Subtract the reagent blank value from the final results or perform a reagent blank adjust.

2- Standard Solution Method

Prepare a 0.25-mg/L trivalent chromium standard as follows:

- 1). Dilute 5.00 mL of Trivalent Chromium Standard Solution, 50-mg/L as Cr^{3+} , to 1000 mL with deionized water. Prepare this solution daily.

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2). To adjust the calibration curve using the reading obtained with the 0.25-mg/L standard solution, touch **Options** on the current program menu. Touch **Standard Adjust**.

3. Touch **On**. Touch **Adjust** to accept the displayed concentration (the value depends on the selected chemical form). If an alternate concentration is used, touch the number in the box to enter the actual concentration, then touch **OK**. Touch **Adjust**.

Method Performance

Precision

Standard: 0.25 mg/L Cr

95% Confidence Limits of Distribution	Program
0.24 – 0.26 mg/L Cr.	100

Standard Additions Method (Sample Spike)

1. After reading test results, leave the sample cell (unspiked sample) in the instrument. Verify the chemical form.

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A summary of the standard additions procedure will appear.

3. Touch **OK** to accept the default values for standard concentration, sample volume, and spike volumes. Touch **Edit** to change these values. After values are accepted, the unspiked sample reading will appear in the top row. See *Standard Additions* in the instrument manual for more information.
4. Snap the neck off a Trivalent Chromium Voluette® Ampule Standard, 12.5-mg/L as Cr³⁺.
5. Prepare three sample spikes. Fill three mixing cylinders (Cat. No. 1896-40) with 25 mL of sample. Use the TenSette® Pipet to add 0.1 mL, 0.2 mL, and 0.3 mL of standard, respectively, to each sample and mix thoroughly.
6. Analyze each sample spike as described in the procedure above, starting with the 0.1 mL sample spike. Accept each standard additions reading by touching **Read**. Each addition should reflect

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approximately 100% recovery.

7. After completing the sequence, touch **Graph** to view the best-fit line through the standard additions data points, accounting for matrix interferences. Touch **Ideal Line** to view the relationship between the sample spikes and the "Ideal Line" of 100% recovery.

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






8. Measurement Procedure (Chromium, Total)

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**Standard Operation Procedure (SOP)
For the Determination of Chromium, Hexavalent**

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Precision

Standard: 0.500 mg/L Cr⁶⁺

95% Confidence Limits of Distribution	Program
0.497 – 0.503 mg / L Cr ⁺⁶	90

3- Standard Additions Method (Sample

Spike)

1. After reading test results, leave the sample cell, or AccuVac Ampul (unspiked sample) in the instrument. Verify the chemical form.
2. Press **OPTIONS>MORE**. Press **STANDARD ADDITIONS**. A summary of the standard additions procedure will appear.
3. Press **OK** to accept the default values for standard concentration, sample volume, and spike volumes. Press **EDIT** to change these values. After values are accepted, the unspiked sample reading will appear in the top row. See the user manual for more information.

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4. Snap the neck off a Chromium Voluette® Ampule Standard, 12.5 mg/L Cr6+.

5. For analysis using powder pillows, use the TenSette® Pipet to add 0.1 mL, 0.2 mL, and 0.3mL of standard, respectively to three 25-mL samples and mix thoroughly. Transfer 10 mL of each solution into a 10-mL sample cell and analyze as described above.

Note: For AccuVac ampules, fill three mixing cylinders* with 50-mL of sample and spike with 0.2 mL, 0.4 mL, and 0.6 mL of standard. Transfer 40 mL from each of the three mixing cylinders to three 50-mL beakers. Analyze each standard addition sample as described in the procedure above.

6. Accept each standard additions reading by pressing **READ**. Each addition should reflect approximately 100% recovery.

7. After completing the sequence, press **GRAPH** to view the best-fit line through the standard additions data points, accounting for matrix interferences. Press **IDEAL LINE** to view relationships between the sample spikes

Chromium Voluette® 12.5 mg/L as Cr⁺⁶
(Ampule Standard)

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.ADJUST / OK

2- Standard Solutions Method

1. To check accuracy, use a 1.0-mg/L Nitrogen Ammonia Standard Solution.
Or, prepare a 1.0-mg/L ammonia nitrogen standard solution by pipetting 1.00 mL of Nitrogen Ammonia Voluette® Ampule Standard, 50-mg/L, into a 50-mL volumetric flask. Dilute to the mark with deionized water. Prepare this solution daily. Perform the Nessler procedure as described above
2. To adjust the calibration curve using the reading obtained with the standard solution, press **OPTIONS>MORE** on the current program menu. Press **STANDARD ADJUSTS**.
3. Press **ON**. Press **ADJUSTS** to accept the displayed concentration. If an alternate Concentration is used, press the number in the box to enter the actual concentration, then press **OK**. Press **ADJUST**

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1. After reading test results, leave the sample cell (unspiked sample) in the instrument.
2. Press **OPTIONS>MORE**. Press **STANDARD ADDITIONS**. A summary of the standard additions procedure will appear.
3. Press **OK** to accept the default values for standard concentration, sample volume, and spike volumes. Press **EDIT** to change these values. After values are accepted, the unspiked sample reading will appear in the top row. See the user manual for more information.
4. Snap the neck off a Nitrogen Ammonia Voluette® Ampule Standard, 50-mg/L NH₃-N.
5. Prepare three sample spikes. Fill three mixing cylinders with 25 mL of sample.

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Use the TenSette® Pipet to add 0.1 mL, 0.2 mL, and 0.3 mL of the 50 mg/L standard, respectively, to each sample and mix thoroughly.

6. Analyze each sample spike as described in the procedure above, starting with the 0.1 ML sample spike. Accept each standard additions reading by pressing **READ**. Each addition should reflect approximately 100% recovery.

7. After completing the sequence, press **GRAPH** to view the best-fit line through the standard additions data points, accounting for the matrix interferences. Press **IDEAL LINE** to view the relationship between the sample spikes and the "Ideal Line" of 100% recovery.

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**Standard Operation Procedure (SOP)
For the Determination of Cyanide Ion
(By Electrode)**

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6. Required Reagents

- 1) Deionized or distilled water for solution and standard preparation.
- 2) Cyanide Ionic Strength Adjuster (ISA): 10M NaOH, (to be prepared from reagent-grade Sodium Hydroxide (NaOH))
- 3) Cyanide Standard: 1,000 mg/L CN⁻, (to be prepared from reagent-grade KCN or NaCN)
- 4) Water deionized

7. Unit of Measurement

Cyanide ions can be measured in units of mg/L (ppm), moles per liter, or any other convenient concentration unit. Table below indicates some concentration unit and conversion factors

mg/L CN ⁻	Moles/L
26.00	1.0×10^{-3}
10.00	3.8×10^{-4}
2.60	1.0×10^{-4}
1.00	3.8×10^{-5}
0.26	1.0×10^{-5}

8. Electrode Storage

The cyanide electrode may be stored for short periods of time in 2.6 mg/L (1.0×10^{-4} M) cyanide solution with ISA added. For longer storage (longer than two weeks), rinse and dry the sensing pellet and cover the membrane tip

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Key procedure:

- ① Pre-preparing sample (dilution if needed)
- ② Cyanide Standard: 1,000 mg/L CN⁻, (to be prepared from reagent-grade KCN .
- ③ Tenfold standard solution ,for example, 20/200,30/300,40/400 starting with higher concentration solution.
- ④ Cyanide Ionic Strength Adjuster (ISA): 10M NaOH, (to be prepared from

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2 – how to get tenfold standard solution:

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**3 - Preparing Ionic standard adjuster ISA
10M (NaOH)**

Fill a 1000 ml beaker with about 900 ml of distilled water. While gently stirring under a hood > slowly add 400 grams of reagent-grade sodium hydroxide.

Transfer the solution quantitatively to a one liter volumetric flask after the solid NaOH has dissolved and beaker has cooled. Dilute to the

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10 M NaOH = 10 mole NaOH

mark with distilled water, cap, upend several times to thoroughly mix the solution.

Store in place a plastic bottle.

TO simplify the preparation, take 250 ml instead of 1000 ml of ISA.

$$\begin{aligned} 10 \text{ M NaOH} &= 10 \text{ mole NaOH} \\ &= 10 \times 40(\text{g- NaOH}) \\ &= 400 \text{ g} \end{aligned}$$

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7-single click on SET UP to chose the parameter to be measured – figure (4)

8 – single click on(Single Channel) icon to enter ion measurement display – figure (5)

Note: use MODE to change between windows.

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
figure 2

Assigning channels for PC6000, PH6500, PC6500 and PCD6500 (example shown is PCD6500)

Figure 3

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- Sample ID:** MANUAL, SEQUENCE, NONE
- Ion Method:** Direct Reading with Standards
- Electrode Type:** AMMONIA
- Measurement Unit:** ppm
- Auto Read Mode:** AUTO, MANUAL
- Ion Stability Criteria:** FAST, MEDIUM, SLOW
- Default Temperature:** 25.0, C, F, K
- Apply Temperature Compensation:** Yes, No; Isopotential Point: 0.0 mV
- Alarm Limits:** ON, Low (10.00 E-7), OFF, High (9.99 E10)
- Print Criteria:** Touch here to edit
- Data Storage Criteria:** Touch here to edit
- Display Criteria:** Touch here to edit

Annotations on the left side of the menu:

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Figure 4


Figure 5

**Standard Operation Procedure (SOP)
For the Determination of Nitrate Ion
(By Electrode)**

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1. Scope and Application: For water, wastewater

2. Summary of Method: Ion Selective Electrode (ISE) (Direct Measurement using ion meter)

3. LIMITS OF DETECTION: the upper limit of detection in pure sodium nitrate solution is 1M. In the presence of other ions, the upper limit of detection is above 10^{-1} M nitrate. But two factors influence this upper limit. Both the possibility of a liquid junction potential development at the reference electrode and the extraction effect influence this limit. Some salts may infuse into the electrode membrane at high salt concentration, causing deviation from the theoretical response. Either dilute samples between 1M and 10^{-1} M or calibrate the electrode at 4or 5 intermediate points. **The lower limit of detection is influenced by the slight water solubility of the ion exchanger used in the sensing portion of the electrode. Nitrate measurements below 10^{-5} M NO_3^- (0.6 ppm as NO_3^-) should employ low level procedures.**

4. Sample Requirements:

All samples must be aqueous and not contain organics which can dissolve in the membrane or extract out the liquid ion exchanger.

The temperature of the sample solutions and of the standard solutions should be the same and below 40°C . About 2% error will be introduced for a 1°C difference in temperature.

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6. Required Reagents

- 1) Deionized or distilled water for solution and standard preparation.
- 2) Nitrate Ionic Strength Adjuster (ISA): 2M $(\text{NH}_4)_2\text{SO}_4$, Code No. (EC-ISA-NT1-BT)
- 3) Nitrate Standard: 1,000 mg/L NO_3^- , (to be prepared from reagent-grade Sodium Nitrate (NaNO_3))
- 4) Water deionized

7 . Unit of Measurement

Nitrate concentrations are measured in units of mg/L (ppm) as Nitrate, moles per liter, or any other convenient concentration unit. Table below indicates some concentration unit and conversion factors

Mg/L NaNO_3	mg/L NO_3^-	Moles/L
850.0	620.0	1.0×10^{-2}
85.0	62.0	1.0×10^{-3}
8.5	6.2	1.0×10^{-4}

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8. Electrode Storage

The nitrate ion electrode may be stored for short periods of time in 620.0 mg/L (1.0×10^{-2} M) nitrate solution. For longer storage (longer than two weeks), rinse and dry the nitrate membrane and cover the tip with any protective cap shipped with electrode. The reference portion of the combination electrode should be drained of filling solution, if refillable and the rubber insert placed over the filling hole.

9. pH Effects

Operating range of the nitrate electrode is from pH2.5 to pH11

Mg/L NaNO ₃	mg/L NO ₃ ⁻	Moles/L
850.0	620.0	1.0×10^{-2}
85.0	62.0	1.0×10^{-3}
8.5	6.2	1.0×10^{-4}

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Procedure

1 - Preparation of Standard Solution 1,000 ppm

(1,000 mg-NO₃⁻/L or ppm) add 1.37 grams of reagent-grade **sodium nitrate (NaNO₃)** to a one liter volumetric flask about half full with distilled water. Swirl the flask gently to dissolve the solid. Full to the mark with distilled water , cap and upend several times to mix the solution .

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To prepare this solution from your own laboratory stock, half fill a 1000 ml volumetric flask with distilled water and add 264 grams of reagent-grade ammonium sulfate, (NH₄)₂SO₄. Swirl the flask gently to dissolve the solid. Fill the flask to the mark with distilled water, cap, and upend several times to mix the contents. ISA is added at the rate of 2 ml of ISA to each 100 ml of standard or sample to adjust the ionic strength to about 0.1 M

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
Figure 3

Assigning channels for PC6000, PH6500,
PC6500 and PCD6500 (example shown is
PCD6500)

The screenshot shows the 'Ion(Ch1) Setup' menu with the following settings and annotations:

- Sample ID:** MANUAL, SEQUENCE, NONE
- Ion Method:** Direct Reading with Standards (Annotation: طريقة الكشف عن الايونات.)
- Electrode Type:** AMMONIA (Annotation: نوع الالكترود)
- Measurement Unit:** ppm (Annotation: وحدة القياس)
- Auto Read Mode:** AUTO, MANUAL (Annotation: ضبط القراءة/ يفضل الالي .)
- Ion Stability Criteria:** FAST, MEDIUM, SLOW (Annotation: معايير استقرار و استجابة الايون.)
- Default Temperature:** 25.0, C, F, K (Annotation: درجة الحرارة.)
- Apply Temperature Compensation:** Yes, No (Annotation: ضبط مقارن لدرجتي حرارة العينة و المحاليل القياسية.)
- Alarm Limits:** ON, Low (10.00 E-7), OFF, High (9.99 E10) (Annotation: حدود الانذار .)
- Print Criteria:** Touch here to edit (Annotation: طباعة المعايير.)
- Data Storage Criteria:** Touch here to edit (Annotation: معايير تخزين المعطيات.)
- Display Criteria:** Touch here to edit (Annotation: عرض المعايير.)

Buttons at the bottom: OK, Cancel, View, Help, Reset.

Figure (4)


Figure (5)

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**Standard Operation Procedure (SOP)
For the Determination of Nitrate
By
Spectrophotometer (DR5000)**

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


4.8 – 5.2 mg/L NO₃⁻-N


3- Standard Additions Method (Sample Spike)

1. After reading test results, leave the sample cell (unspiked sample) in the instrument. Verify the chemical form.
2. Press **OPTIONS>MORE**. Press **STANDARD ADDITIONS**. A summary of the standard additions procedure will appear.
3. Press **OK** to accept the default values for standard concentration, sample volume, and spike volumes. Press **EDIT** to change these values. After values are accepted, the unspiked sample reading will appear in the top row.
4. Open a bottle of Nitrate Nitrogen Standard, 100-mg/L NO₃⁻-N.
5. Prepare three sample spikes. Fill three sample cells with 10 mL of sample. Use the micropipette to add 0.1 mL, 0.2 mL, and 0.3 mL of standard, respectively, to each sample and mix thoroughly.
6. Analyze each sample spike as described in the procedure below, starting with the 0.1 mL sample spike. Accept each standard additions reading by pressing **READ**. Each addition should reflect approximately 100% recovery.
7. After completing the sequence, press **GRAPH** to view the best-fit line through the standard additions data points, accounting for matrix interferences. Press **IDEAL LINE** to view Relationships between the sample spikes and the “Ideal Line” of 100% recovery.

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**Standard Operation Procedure (SOP)
For the Determination of Nitrite
By
Spectrophotometer (DR5000)**

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


Date: Aug. 2007

Page 3 of 5

<p>(1) - Reagent blank</p> <p>For more accurate results, determine a reagent blank value for each new lot of reagent. Follow the procedure using deionized water instead of the sample. Subtract the reagent blank value from the final results or perform a reagent blank adjust.</p> <p>Method Performance</p> <p>Precision</p> <p>Standard: 0.150 mg/L NO₂⁻-N</p> <table border="1" data-bbox="167 862 829 990"><tr><td data-bbox="167 862 829 936">95% Confidence Limits of Distribution</td></tr><tr><td data-bbox="167 936 829 990">0.147 – 0.153 mg/L NO₂⁻-N</td></tr></table>	95% Confidence Limits of Distribution	0.147 – 0.153 mg/L NO₂⁻-N	<p>1- طريقة معايرة الشاهد :</p>
95% Confidence Limits of Distribution			
0.147 – 0.153 mg/L NO₂⁻-N			

6. Measurement Procedure (Phosphorus, Reactive)

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Step	Operation	Remarks (ملاحظات)	التشغيل
7	Wipe the blank and insert it into the cell holder with the fill line facing the user. Press ZERO . The display will show: 0.000 mg/L NO ₂ -N		
8	Wipe the prepared sample and insert it into the cell holder with the fill line facing the user. Results are in mg/L NO ₂ -N		

**Standard Operation Procedure (SOP)
For the Determination of Chloride Ion
(By Electrode)**

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5. Sample Requirements:

All samples must be aqueous and not contain organics which can dissolve the epoxy electrode body and/or the cement bonding the sensing crystal to the electrode body. Inorganic solutions will not affect the electrode. Infrequent measurements in solutions containing methanol, ethanol, benzene, and acetonitrile are permitted. Highly polar solvents slowly attack the electrode.

The temperature of the standard solutions and of the sample solutions should be the same and below 50°C.

The pH range for the chloride ion electrode is 2 – 12. Neutralize samples outside this range with acid or base to bring them in range.

6. Necessary Equipment and Supplies:

- 1) Ion meter (CyberScan PCD6500 Meter)
Eutech Chloride Ion Combination
Epoxy-body Electrode, Code no.
EC-CLO-03.
- 2) Volumetric Flask (for preparation of standard solutions and ISA)
- 3) Beaker (150 – 250 mL)
- 4) Magnetic stirrer
- 5) Polishing Paper: to polish dirty or etched electrode membranes.

7. Required Reagents

- 1) Deionized or distilled water for solution and standard preparation.
- 2) Ionic Strength Adjuster (ISA)

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- 3) Chloride Standard: 1,000 mg/L F⁻, (to be prepared from reagent-grade sodium chloride NaCl)
- 4) Water deionized

8. Unit of Measurement

Chloride concentrations are measured in units of mg/L (ppm), moles per liter, or any other convenient concentration unit. Table below indicates some concentration unit and conversion factors.

mg/L Cl ⁻	Moles/L
354.50	1.0×10^{-2}
35.45	1.0×10^{-3}
3.55	1.0×10^{-4}

9. Electrode Storage

The electrode may be stored for short periods of time in 1.0×10^{-2} M chloride solution. For longer storage (longer than two weeks), rinse and dry the sensing pellet and cover the membrane tip with any protective cap shipped with the electrode. The reference portion of the combination electrode (or the outer chamber of the reference electrode) should be drained of filling solution, if refillable, and the rubber insert placed over the filling hole.

10. pH Effects

Hydroxide ion interferes with measurements of low levels of chloride although the electrode can be used over a reasonable pH range >

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Key steps to conduct the Procedure

- ① Pre-preparation of sample(preparing dilutions if needed , for instance , dying factory.)
- ② Preparation of Standard Solution 1,000 ppm (1,000 mg-Cl⁻/L) using reagent-grade **sodium chloride (NaCl)**.
- ③ Prepare standard chloride solutions whose concentrations vary by tenfold (20/200,30/300,40/400) and high concentration standard solution should be prepared firstly.
- ④ Preparing Ionic standard adjuster ISA (5MNaNO₃).
- ⑤ Electrode Slope Check (using standard chloride solutions whose concentrations vary by tenfold and ISA (5MNaNO₃))
- ⑥ Direct measurement.

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11. Preparation of Standard Solution and ISA

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Preparation of 1,000 ppm (1,000 mg-Cl⁻/L):

Add 1.65 grams of reagent-grade **sodium chloride (NaCl)** to a one liter volumetric flask about half-full of distilled water. Swirl the flask to dissolve the solid. Fill the flask to the mark with distilled water, cap, and upend the flask several times to mix the contents.

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I – switch on power supply , then , connect the target ion electrode directly to Channel 1 or 2 (to

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

figure (2)

Assigning channels for PC6000, PH6500,PC6500 and PCD6500 (example shown isPCD6500)

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
- Sample ID:** MANUAL, SEQUENCE, NONE
- Ion Method:** Direct Reading with Standards
- Electrode Type:** AMMONIA
- Measurement Unit:** ppm
- Auto Read Mode:** AUTO, **MANUAL**
- Ion Stability Criteria:** FAST, **MEDIUM**, SLOW
- Default Temperature:** 25.0, C, F, K
- Apply Temperature Compensation:** Yes, No, Isopotential Point: 0.0 mV
- Alarm Limits:** ON, Low: 10.00 E-7, OFF, High: 9.99 E10
- Print Criteria:** Touch here to edit
- Data Storage Criteria:** Touch here to edit
- Display Criteria:** Touch here to edit

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Figure (5)

**Standard Operation Procedure (SOP)
For the Determination of Fluoride Ion
(By Electrode)**

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1. **Scope and Application:** For water, wastewater

2. **Summary of Method:** Ion Selective Electrode (ISE) (Direct Measurement using ion meter)

3. **Measurement range (Limits of Detection):**

Fluoride concentration down to $1.0 \times 10^{-6} \text{M}$ (0.02 mg/L) fluoride can be measured in neutral solutions. Since sample contamination can be a factor in low level fluoride measurements, care must be taken in making determinations below $1.0 \times 10^{-5} \text{M}$. The upper limit of detection is a saturated fluoride solution.

4. **Sampling and Storage**

Preferably use polyethylene bottles for collecting and storing samples for fluoride analysis. Glass bottles are satisfactory if previously they have not contained high-fluoride solutions. Always rinse bottle with a portion of sample.

5. **Sample Requirements:**

All samples must be aqueous and not contain organics which can dissolve the epoxy electrode body and/or the cement bonding the sensing crystal to the electrode body. Inorganic solutions will not affect the electrode. Infrequent measurements in solutions containing methanol, acetone, or dioxin are permitted. Highly polar solvents, such as CHCl_3 or DMF, should not be contained in the samples.

The addition of Total **Ionic Strength Adjuster Buffer** (TISAB) to samples and standards will

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adjust the pH to 5.0-5.5. Samples must be above pH 5 to avoid forming complexes with hydrogen ions and below pH 7 to avoid interference by hydroxide ions.

The temperature of the standard solutions and of the sample solutions should be the same and below 80°C. The use of TISAB 1 (for high level) also preferentially forms complexes with aluminum and with iron, breaking the complexes that fluoride forms with these ions. With 1 ppm fluoride present, up to 3-5 ppm aluminum or iron is complexed. If higher levels of aluminum or iron are present, use TISAB3.

6. Necessary Equipment and Supplies:

- 1) Ion meter (Cyber Scan PCD6500 Meter)
- 2) Fluoride Ion Selective Electrode (Eutech Cyanide Ion Combination Epoxy-body Electrode, Code no. EC-FO-03)
- 3) Volumetric Flask (for preparation of standard solutions and ISA)
- 4) Beaker (150 – 250 mL)
- 5) Magnetic stirrer
- 6) Polishing Paper: to polish dirty or etched electrode membranes.

7. Required Reagents

- 1) Deionized or distilled water for solution and standard preparation.
- 2) Total Ionic Strength Adjuster Buffer (TISAB 1 and TISAB 2)

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
11. Electrode Storage

The fluoride electrode may be stored for short periods of time in 190 mg/L (1.0×10^{-2} M) fluoride solution with TISAB added. For longer storage (longer than two weeks), rinse and dry the sensing pellet and cover the membrane tip with any protective cap shipped with electrode. The reference portion of the combination electrode (or the outer chamber of the reference electrode) should be drained of filling solution, if refillable, and the rubber insert placed over the filling hole. The fluoride electrode should never be stored in distilled water.

12. pH Effects

Hydrogen complexes a portion of fluoride in solution, forming the un-dissociated acid HF and the ion HF_2^{-1} in acid solutions with a pH below 5. The proportion of free fluoride ion in acid solutions is shown in Figure below.

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
Fraction of free F^- as a function of solution pH, where hydrogen is the only complexing species

When the level of hydroxide is greater than one-tenth the level of fluoride ion present, hydroxide ion interferes with electrode response to fluoride. As an example, no hydroxide interference with fluoride measurements take place at pH 7 when the hydroxide concentration is 1.0×10^{-7} or less. As the pH increases, the hydroxide interference becomes appreciable. At pH 10, hydroxide ion concentration is $1.0 \times 10^{-4} M$ and no error is found in measurements of $1.0 \times 10^{-2} M$ fluoride. At the same hydroxide ion concentration and a fluoride concentration of $1.0 \times 10^{-4} M$, about a 10% measurement error appears. At a fluoride concentration of $1.0 \times 10^{-5} M$, considerable error exists in a pH 10 solution. Following Figure illustrates these errors.


The addition of TISAB 1 or TISAB 2 to all fluoride samples and standards buffers the pH between 5.0-5.5 to help avoid hydroxide interference or the formation of hydrogen complexes of fluoride

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enter ion measurement display – figure (5)

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Ion(Ch1) Setup pH(Ch2) Setup pH(Fet) Setup Cond Setup

Sample ID MANUAL SEQUENCE NONE

Ion Method Direct Reading with Standards

Electrode Type AMMONIA

Measurement Unit ppm

Auto Read Mode AUTO MANUAL

Ion Stability Criteria FAST MEDIUM SLOW

Default Temperature 25.0 C F K

Apply Temperature Compensation Yes No Isopotential Point 0.0 mV

Alarm Limits ON Low 10.00 E-7 OFF High 9.99 E10

Print Criteria Touch here to edit

Data Storage Criteria Touch here to edit

Display Criteria Touch here to edit

OK Cancel View Help Reset

Figure 4


Figure 5

**Standard Operation Procedure (SOP)
For the Determination of Sulfide**

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sulfate to form methylene blue. The intensity of the blue color is proportional to the sulfide concentration . high sulfide levels in oil field waters may determined after proper dilution. Test results are measured at 665 nm.

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Collect sample in clean plastic or glass bottles. Fill completely and cap tightly. Avoid excessive agitation or prolonged exposure to air. Analyze samples immediately.

Don't store sample for more than 7 days before analysis, and pH to be set at least 9.

For sample storage, follow below:

Add 2N Zinc acetate 4 drops before sample collection for every 100mL sample. Add then adjust pH by using Sodium Hydroxide to be more than 9.

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Preparation of Sulfide Standards

Take care in preparing reliable stock solutions of sulfide for calibration and quality control. Prepare sulfide standards from sodium nonahydrate ($Na_2S \cdot 9H_2O$) crystals.

Preferably remove single crystals of $Na_2S \cdot 9H_2O$ from reagent bottle with nonmetallic tweezers; quickly rinse in degassed reagent water to remove surface contamination. Blot crystal dry with a tissue, then rapidly transfer to a tared, stoppered weighing bottle containing 5

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to 10 mL degassed reagent water. Repeat procedure until desired amount of sodium sulfide is in weighing bottle. Determine amount of $\text{Na}_2\text{S} \cdot 9\text{H}_2\text{O}$ in weighing bottle by difference, then multiply the weight by 0.133 to determine the amount of S^{2-} . Avoid excess agitation and mixing of the solution with atmospheric oxygen. Quantitatively transfer and dilute entire contents of weighing bottle to an appropriate size volumetric flask with degassed reagent water to prepare a known concentration sulfide stock solution (3.750 g $\text{Na}_2\text{S} \cdot 9\text{H}_2\text{O}$ diluted to a final volume of 500 mL will give a stock solution of which 1.00 mL = 1.00 mg S^{2-}) Store stock solution with minimum headspace for no more than 1 week.

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



**Standard Operation Procedure (SOP)
For the Determination of total suspended solids**

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
necessary to increase the diameter of the filter or decrease the sample volume.

a. Preparation of glass-fiber filter disk

- ① Place a filter disk on the filter holder. Assemble filter holder in suction flask apparatus, connect to vacuum source and apply vacuum.
- ② Wash the filter disk with three successive 20-mL portions of distilled water. Continue to apply vacuum for 2-3 minutes after the water has passed through the filter. Discard the filtrate.
- ③ Remove the filter from filtration apparatus and transfer to an inert aluminum weighing dish. Dry in an oven at 103 – 105°C for at least 1 h.
- ④ Cool in desiccator to balance temperature and weigh on an analytical balance.
- ⑤ Repeat the cycle of drying, desiccating and weighing until a constant weight is obtained or until weight change is less than 4% of the previous weighing or 0.5 mg, whichever is less.
- ⑥ Store in desiccator until needed.

b. Sample analysis

- ① Remove the filter disk from the Desiccator, weigh it and record its weight.
- ② Place the filter in filter holder and assemble the filter holder in the suction flask apparatus. Connect to the vacuum source and apply vacuum.
- ③ Wet the filter with a few drops of distilled water to seat the filter.
- ④ Shake the sample vigorously and

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

Procedure of Determination of SS (Gravimetric method)

Short Title: Total Coliform (Membrane Filter Technique)

Revision No.: 0

Date: Sept. 2007

Page 1 of 27

**Standard Operation Procedure (SOP)
For the Determination of
Total Coliform (Membrane Filter Technique)**

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What is Total Coliform? (1)

- ❁ The term coliform bacteria represents a vaguely defined group of organisms which have a long history in water quality assessment.

What is Total Coliform? (2)

- ❁ Pathogenic organism \neq Coliform bacteria
- ❁ Number in water
 - Pathogenic organism: Small
 - Coliform bacteria: Large
- ❁ Coliform bacteria mostly includes pathogenic organism
- ❁ Test (Measurement)
 - Pathogenic organism: Difficult
 - Coliform bacteria: Relatively easy

What is Total Coliform? (3)

- The term “total coliform” refers to **a large group** of:
 - Gram-negative,
 - Rod-shaped bacteriathat **share several characteristics**.


**Pathogens: bacteria, protozoa, and viruses
that make people sick**

What is Total Coliform? (4)

- Non-disease-causing organisms found in soil or vegetation and in the intestinal tract of warm-blooded animals (fecal coli.).
- Present in much larger numbers than the more dangerous pathogens, and react to the natural environment and treatment processes in a manner and degree similar to pathogens.

What is Total Coliform? (5)

- One of the general sanitary water quality indicator that suggests the possibility of presence of pathogenic organism.

Why test for coliform bacteria? (1)

- Most coliforms are not pathogens
- The presence of very few coliforms in water
➡ Water probably contains no pathogenic organisms
- Presence of large number of coliforms
➡ Very high probability of contamination by pathogenic organism
- Number of pathogenic organism ➡ Relatively small
➡ Very difficult to isolate and identify specific pathogenic organism

Why test for coliform bacteria (2)

5. Observing/Testing of coliform bacteria

- ➔ Increase or decrease of many pathogenic organism can be estimated
- ➔ Total coliforms are indicators and are more common and easy to test

Total coliforms are mostly natural residents of soil and water. Coliform bacteria are those that are usually found in the fecal material of animals. Their presence usually means that the water may be contaminated by sewage effluent. Finding the source of the problem and correcting it is very important.

Where they are found? (1)

- There are many pollution sources.
- Domestic animals contribute heavily to the population of coliform
- Including runoff from:
 - Woodland
 - Pastures
 - Feedlots
 - Septic tank
 - Sewage plants
 - Animals and wild fowl
 - Others

How to measure/test coliforms?

- ◆ In the laboratory, tests are to conduct using following principles
 - Grown in or on a medium containing lactose, at a temperature 35 or 37°C
 - Provisionally identified by the production of acid and gas from the fermentation of lactose

Selection of Analytical Method

- Commonly used two techniques
 1. “Multiple fermentation tube” technique or “Most probable number” technique

Measured portion of water sample are placed in test-tubes containing a culture-medium. The tubes are then incubated for a standard time at a standard temperature


2. “Membrane filter” technique ← **To be adopted**

Measured volume of sample is passed through a fine filter that retains bacteria. The filter is then placed on culture medium and incubated.

Comparison of Method

Multiple fermentation tube technique	Membrane Filter technique
Slower: Requires 48 hrs for a positive	More rapid: Requires 24 hrs
More labor- intensive	Less labor-intensive
Requires more culture medium	Requires less culture medium
More sensitive	Less sensitive
Result obtained indirectly by statistical approximation (low precision)	Results obtained directly by colony count (high precision)
Not readily adaptable for use in the field	Readily adapted for use in the field
Applicable to all types of water	Not applicable to turbid water

Size of Variety of Particles in Water

Membrane Filter Technique

- ❁ Introduced in the late 1950s as an alternative to the “Multiple fermentation tube” technique.
- ❁ Offers the advantage of isolating discrete colonies of bacteria.

- ❁ Whereas the multiple fermentation tube technique only indicate the presence or absence of an approximate number or organism.

Procedure of Membrane Filter Technique (Outline)

- A definite volume of sample; in the case of drinking water normally 100 mL, is passed through a **47 mm** membrane of uniform pore diameter, usually **0.45 μ** , using a filter funnel and vacuum system.
- Any organism in the sample are **trapped/retained on the surface of the membrane**. The filter is then placed in a petri dish with **nutrient medium** and incubate at an appropriate temperature.
- The passage of nutrients through the filter facilitates the growth of organisms on the upper surface of the membrane.
- The discrete colonies that form on the surface of the membrane is transferred to a colony counter and number of colonies is to be counted.

Outline Procedure


- ① Measured volume of water is filtered through a cellulose acetate membrane of uniform pore diameter
- ↓
- ② Bacteria are retained on the surface of the membrane
- ↓
- ③ The membrane is placed on a suitable selective medium (culture media) in a sterilized container, and incubated at an appropriate temperature
- ↓
- ④ If coliforms and/or faecal coliforms are present in the water sample, characteristic colonies form that can be counted directly

(Note: All materials and equipment must be sterilized prior to use)

Apparatus and Materials

- ① Dilution bottles or tubes,
- ② Pipets and graduated cylinders,
- ③ Containers for culture media,
- ④ Culture dishes (Petri-type dishes),
- ⑤ Filtration units,
- ⑥ Filter membrane,
- ⑦ Forceps,
- ⑧ Colony counter,
- ⑨ Sterilizing device (autoclave, oven, etc.)

Procedure of Membrane Filter Technique


Colony Counting (1)

- ⚙ When counting the colonies the entire surface of the filter should be scanned using a 10x -- 15x binocular, wide-field dissecting microscope, etc.
- ⚙ Number of colonies yield
 - Ideally About 50/membrance
 - Not more than 200/membrance
 - Preferably 20 to 80/membrance
- ⚙ Counts for each filter should be recorded on the laboratory data sheet.

Colony Counting (2)

Colony Counting
Technique

Culture Media (1)

- A substance or material used for the growth of coliform bacteria
- To be selected in accordance with kind of coliform bacteria detected
- Can be obtained in the form of dehydrated (Granular culture media)
- M-Endo is a representative culture media for the enumeration of Escherichia coil (e-coli) in the Standard Total Coliform Membrane Filter Method

Culture Media (2)

- ❁ Heat-sensitive
- ❁ Don't heat any longer than necessary
- ❁ Storage of dehydrated culture media
 - To be stored in a dry, dark place at a temperature of about +15 to 25°C
 - Containers should be well sealed and tightly closed
 - Absorption of water leads to pH shifts and eventually clumping

Selection of Sample Size

- ❁ To be governed by expected bacteria density
- ❁ Ideal sample volume
 - About 50 coliform colonies per membrane,
 - Not more than 200 colonies,
 - Preferably 20 – 80 colonies per membrane
- ❁ To be selected in accordance with sample type

Suggested volume to be filtered (1)

Sample type	Sample volume (mL)					
	100	10	1 ¹	0.1 ^{1,2}	0.01 ^{1,2}	0.001 ^{1,2}
Treated drinking water	★					
Partially treated drinking water	★	★				
Recreational water			★	★		
Protected source water		★	★			
Surface water			★	★		
Wastewater			★	★	★	
Discharge from sewage treatment plant		★	★	★	★	
Ponds, rivers, stormwater runoff			★	★	★	★
Raw sewage				★	★	★
Feedlot runoff				★	★	★
Well, spring	★	★				

Suggested volume to be filtered (2)

- ¹ Small volume should be added to the filtration apparatus together with a minimum of 9 mL of sterile diluent to ensure adequate dispersal across the surface of the filter membrane.
- ² 1.0, 0.1, 0.01 and 0.001 mL volumes are filtered after first preparing serial dilutions of the sample.

To filter:

- 1.0 mL of sample, use 10 mL of 1/10 dilution
- 0.1 mL of sample, use 10 mL of 1/100 dilution
- 0.01 mL of sample, use 10 mL of 1/1,000 dilution
- 0.001 mL of sample, use 10 mL of 1/10,000 dilution

Sample Collection, Preservation, and Handling

- ❁ Clean all glassware thoroughly with a suitable detergent and hot water, rinse with hot water to remove all traces of residual washing compound, and finally, rinse with distilled water.
 - ❁ Sterilize glassware for not less than 2 hr at a temperature of 170°C.
 - ❁ Sterilize sample bottles not made of plastic, as above, or in an autoclave at 121°C for 15 min.
- 

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Standard Operation Procedure (SOP)
For the Determination
Of
Chemical Oxygen Demand (COD)
By
Open Reflux Method

Prepared by: _____ **Date:** _____
Chemist

Reviewed by: _____ **Date:** _____
Laboratory Chief or Quality Assurance Coordinator

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
<p>1. Scope and Application: For water, wastewater</p> <p>2. Summary of Method: Open Reflux Method (Oxidizer: Potassium dichromate, $K_2Cr_2O_7$)</p> <p>3. Measurement range: $>50 \text{ mgO}_2/\text{L}$ and $< 50 \text{ mgO}_2/\text{L}$</p> <p>4. Equipment Required:</p> <ol style="list-style-type: none">① 500-milliliter (ml) or 250-milliliter (ml) Erlenmeyer flask with standard tapered glass joints② Liebig reflux condensers (12-inch) with standard tapered glass joints③ Electric hot plate or mantle heater④ Volumetric pipets (10, 25, and 50-ml capacity)⑤ Buret, 50 ml or 25 ml - 0.1 ml accuracy⑥ Buret stand and clamp⑦ Analytical balance, accuracy 0.001 gram (g)⑧ Spatula⑨ Volumetric flasks (1,000 ml capacity)⑩ Boiling beads, glass⑪ Magnetic stirrer and stirring bars or glass beads <p>5. Reagents/Chemicals Required</p> <ol style="list-style-type: none">① Potassium dichromate ($K_2Cr_2O_7$) 0.25 N or 0.025N② Sulfuric acid (H_2SO_4, $d = 1.84$) silver sulfate (Ag_2SO_4) solution (Sulfuric acid reagent)③ Mercuric sulfate ($HgSO_4$) crystals④ Ferrous ammonium sulfate (FAS) $[Fe(NH_4)_2(SO_4)_2] \cdot 6H_2O$, 0.25N or	<p>1. المجال والتطبيق : <u>الماء أو 1:1 حمص الكبريت</u> <u>لضبط الـ</u></p>
---	---

0.025N

- ⑤ Ferriin indicator (1, 10-phenanthroline and ferrous ammonium sulfate)
- ⑥ Potassium hydrogen phthalate (KHP)
- ⑦ Conc. Sulfuric acid
- ⑧ Distilled water

6. Principle of Determination

- Most types of organic matter are oxidized by a boiling mixture of chromic and sulfuric acids,
- Sample is refluxed in strongly acid solution (digestion) with a known excess of potassium dichromate ($K_2Cr_2O_7$),
- After digestion, the remaining unreduced $K_2Cr_2O_7$ is titrated with ferrous ammonium sulfate to determine the amount of $K_2Cr_2O_7$ consumed and the oxidizable organic matter is calculated in terms of oxygen equivalent.

6. Sample Collection, Storage, and Preservation

Test unstable samples without delay. If delay before analysis is unavoidable, preserve sample by acidification to $pH \leq 2$ using conc. H_2SO_4 .

Preferably acidify any sample that cannot be analyzed the same day it is collected. Blend samples containing settleable solids with a homogenizer to permit representative sampling.

Make preliminary dilutions for wastes

containing a high COD to reduce the error inherent in measuring small sample volume.


Refrigerate

➤ **Preparation of Reagents/Chemicals**

1. Potassium dichromate (K₂Cr₂O₇) [Standard potassium dichromate: 0.25 N]

	Operation	Remarks	Operation (Arabic)
P-1	Weigh about 15g of Potassium dichromate (K ₂ Cr ₂ O ₇) and dry in oven at 103°C for 24 hours.		
P-2	After drying, dissolve 12.259g Potassium dichromate (K ₂ Cr ₂ O ₇) in distilled water to 1.0 liter volume in a volumetric flask.		
Note			
M.W. of Potassium dichromate = $39.1 \times 2 + 52.0 \times 2 + 16 \times 7 = 294.2$ This reagent undergoes a six-electron reduction reaction; the equivalent concentration is $6 \times$ $12.259 \div 294.2 = 0.04167 \text{ M}$ $0.04167 \times 6 = 0.25 \text{ N}$ (1 mL solution is equivalent to 2 mg O ₂)			

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Recording Format for Measuring COD_{Cr} by Open Reflux Method

Sample Information

No.	Governorate	Sample Code	Sampling Date	Type/Kind of Sample	Sample Pre-Treatment	Note
1						
2						
3						
4						
5						

Standardization of FAS

1	Standard Pottasium dichromate (V _{PD}):	10.0	mL
2	Addition of conc. H ₂ SO ₄ :		mL

Titration:

Graduation of Burette (Start):		mL
Graduation of Burette (End):		mL
Amount of Titrant (V _{FAS}):	10.00	mL

Calculation of Normality:

$$\text{Normality of FAS solution} = V_{PD}/V_{FAS} \times N_{PD}$$

Normality: _____ N

Reagent/Chemicals prepared

1	Standard Potassium dichromate (K ₂ Cr ₂ O ₇) (N _{PD}):	0.250	N
2	Ferrous ammonium sulfate (FAS) (Fe(NH ₄) ₂ (SO ₄) ₂ ·2H ₂ O):	0.250	M/N

Measurement of COD_{Cr} by Open Reflux Method

Reflux No.	Sample Code/No.	Volume of Sample taken to Reflux	Dilution	Treatment	Standard Potassium dichromate	Mercuric Sulfate HgSO ₄	Suifuric acid reagent
A		mL			mL	g	mL
B		mL			mL	g	mL
C		mL			mL	g	mL
D		mL			mL	g	mL

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Titration of Excess K₂Cr₂O₇ with FAS

Reflux No.	Sample Code/No.	Amount of FAS titrated		Graduation of Burette		Result of Calculation
				Start	End	
A		mL	←	mL	mL	mgO ₂ /L
B		mL	←	mL	mL	mgO ₂ /L
C		mL	←	mL	mL	mgO ₂ /L
D		mL	←	mL	mL	mgO ₂ /L

Calculation of COD

$$\text{COD as mg O}_2/\text{L} = [(A - B) \times M \times 8000] / (\text{mL Sample})$$

where:

A = mL of FAS used for blank,

B = mL of FAS used for sample,

M = Normality of FAS

8,000 = milliequivalent weight of oxygen × 1,000mL/L

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Standard Operation Procedure (SOP) For Water purifying (De-Ionizing)

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◆ Refer the "SOP for Water Stills"
at the same time

1. Scope and Application:

For Generating De-ionized water that is necessary for Atomic Absorption Spectro-Photometer and other experiments where de-ionized water use is necessary

2. Summary of Method:

Water purification system using ion-exchange plus distillation.

3. Raw Water Requirements:

In order to extend the lifetime of ion-exchanger, raw water hardness should not be very high, therefore, distilled water use is advised as the raw water of the de-ioniser.

The temperature of raw water should be below 40C. As the high temperature may reduce the lifetime of the ion-exchanger, and will critically affect the performance of water condenser.

4. Necessary Equipment and Supplies:

- 1) Raw (distilled) water barrel
- 2) Container for submerged pump
- 3) Submerged pump
- 4) Ion-Exchanger
- 5) EC meter
- 6) Distiller
- 7) De-ionized water Tank
- 8) Excess Water Tank
- 9) Set of rubber pipes, clamps

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
Figure 1: Composition of Water Purifier

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6. Safety Measures

1- During the process, high temperature is generated by the condenser, hot water is discharged to de-ionized water container, and excess water container, therefore, it is advisable to wear thermally isolated-waterproof gloves if contact to hot apparatus or containers is necessary.

2- avoid switching electric current with wet hands (or gloves)

The figure below illustrates hot points of the apparatus

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6 . Countermeasures in case of burns due to contact to hot substances

1- in case of severe hand burns

- gently remove clothes (gloves, roll up trousers or shirt... etc.)
- immediately rinse with cold tap water for at least 2 minutes , to allow the cooling of the deep skin layers (dermis).
- if possible apply ice (avoid direct contact to ice with a piece of clean cloth such as a towel), this will locally anesthetize the burns.
- don't use any antiseptics
- don't place any bandages
- apply proper ointment for burns if available
- seek medical advise.

7. Countermeasures in case of severe electric Shock

- immediately switch off electricity from the main panel, if you didn't know how, ask some one who knows!
- Check if you are alone. If there are other people around, instruct them to call an ambulance right away. Or ask them to assist you to deliver the injured person to the hospital quickly, and proceed doing so after you confirm that the area is secured (power is cut);
- Check for a response and breathing, yell to the person to see if he is conscious. At this stage, do not touch the victim. If you are satisfied that it is safe, start

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resuscitating the victim. If you have no first aid knowledge, skip to the next step if someone has not already done it for you.

4. Call **emergency services** for an ambulance.
5. If the breathing and pulse are steady, attend to injuries. Cool the burns and cover with dressings that won't stick. Never put ointments or oils onto burns.
6. Talk calmly and reassuringly to the conscious victim until the ambulance arrives.

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Concept of the Process

Figure 3 illustrates water flow directions through the De-Ionizing cycle, where as figure 4 is the schematic diagram

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- 1- Raw water flow by gravity to the submerged pump tank, then pumped through the ion exchange unit and then to the distiller
- 2- At the distiller water is firstly used for cooling the condenser, then falls to the boiler where it gets evaporated, and condensed water is collected and flow by gravity to the De-Ionized water tank
- 3- Excess water of the cooling-boiling process flow by gravity to the excess water tank
- 4- Figure 5 illustrates the concept of De-Ionizing Distiller

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Concept Of The De-Ionizing Distiller

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- 1- After passing the ion exchange unit, water continue to flow by pressure to the condenser cooling water inlet (7)
- 2- After serving as cooling water, water is collected at the buffer tank (6)
- 3- When water pressure increase at the buffer tank, the boiler inlet valve opens, and allows pressure to be equalised and

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
to keep the coil (3) submerged.

- 4- Water evaporates and condenses, then collected at the De-Ionized water outlet (10)
- 5- When water level at the boiler (10) exceeds its capacity, excess water discharge valve (7) opens and discharge all excess water from the buffer tank and consequently; the boiler.
- 6- If for, any reason, excess water valve fail to open, excess water will be discharged from the overflow (9)

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Procedure

1 – Confirm Water Level in the Raw Water

Tank

The more water you have here, the more water you may generate, as only about 1/10th of raw water is converted to De-Ionized, and the remaining is recycled as excess water; i.e. raw water quantity in this example is barely enough to generate 1 Litre

2- Confirm raw water Temperature

Water temperature should not exceed 40°C, as it may damage the ion-exchange unit, the submerged pump tank, and will drastically reduce the condenser performance and thus will reduce the productivity below feasibility point.

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switch, and each requires considerable capacity, you cannot use the simultaneously.

Hint: *you may save some time and raw water by moving forward step 8 to follow step 2, if only you can perform 3, 4, 5, 6, 7, 8 and 9 before the evaporation of the remaining water at the boiler.*

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10- regulate the water flow;

a. limit the flow through the control valve;

boiler's productivity is limited, so you can't increase productivity by supplying more water than what is actually needed.

- (1) Close the valve, and slightly open it, if the water level at the boiler decreases below the minimum point and continue to decrease, open it slightly, and wait for 1 minute
- (2) If the water level starts to recover, and remained within the green range, then, you have reached the optimum point
- (3) If the water level quickly increases and excess water is frequently discharged, that means that you may limit the flow a little bit more.
- (4) If water level remained within the green range, and water is not frequently discharged through the excess water gate, then the flow rate is acceptable, and you may safely leave the unit unwatched for some time.

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
If excess water pipe curves for any reason (like thermal expansion) it will cause the malfunction of the excess water discharge valve, and water may frequently overflow in the boiler. Keep the pipe as straight as possible.

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b. adjust the flow of the raw water supply valve

by monitoring the tempo of the increase/decrease of water level in the submerged pump tank, you may adjust the flow of the raw water to equalise the pumping speed.

c. comments

keep watching the water level at the boiler, submerged pump tank, and raw water tank, and many factors may affect your adjustments such as; change of water flow from supply tank (due to level decrease), changes of boiler performance due to electrical instability... etc

d. Countermeasures

- in case of power failure, proceed shut down

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tank to be filled with distilled water, the operation of De-Ionizing kit could not be escaped using any available amount, in this case, water may be allowed to be circulated during the process, but temperature of raw water should be monitored, and circulation should be stopped if the temperature of raw water exceed 40°C

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**Standard Operation Procedure (SOP)
For Calibration of Hach DR5000 Spectro-photometer
by creating User Programs using the method of
reading known concentration of Standard Solutions at a defined wave length
1- General Manual of the Method**

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Refer the "SOP for Operation of DR5000 Spectrophotometer" and SOPs of Determination the parameters using DR5000

1. Scope and Application:

To carry out "made to measure" analyses. This is useful when the preparation of standard solutions of concentrations within the measurement limit of the available reagents is possible; i.e measurements are reliable.

To secure all factors mentioned above, standard solutions should be of good quality; i.e valid, or prepared from good quality chemicals, reagents are valid and reliable, glassware and tools used for preparation of standard solution are clean and precise... etc.

The method is particularly important when it is not possible to procure Hach reagents, or compatible. And when it is deemed to be necessary to re-calibrate the device for a different kind of reagents or when a different method is applied rather than Hach standard ones.

2. Summary of Method:

Calibration (derivation of Concentration/ Absorption Graph) of Hach DR5000 spectro-photometer by reading a series of known concentration standard solutions..

3. Necessary Equipment and Supplies:

- 1) DR 5000 Spectrophotometer
- 2) Standard Cells (10ml, 25 ml... etc.)
- 3) Glassware; pipettes, micropipettes, beakers, volumetric flasks, cylinders... volumes and precision is judged on the basis of the desired parameter type, range of measurement, and method.
- 4) Miscellaneous apparatus according to the

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<p>method of preparation and testing of the standard solution; e.g. heater, water bath... etc.</p> <p>5) Calculator</p> <p>6) Registration forms</p>	
<p>4. Necessary Chemicals and Reagent</p> <p>1- Distilled water, but preferable De-ionised water, as the process is quite sensitive.</p> <p>2- Suitable reagents for the parameter / method</p> <p>3- Standard Solutions of suitable concentration (within the range of measurement)</p>	<p>4. إجراءات السلامة</p>
<p>5. Principle of the process:</p> <p>The principle of the process is the core of the concept of the spectrophotometer; i.e. linking between absorbance of the tested sample at a pre-determined wave length (visible or UV) and its concentration.</p> <p>To conduct this calculation, the preinstalled programmes of the spectrophotometer conduct correlation via interpolation (usually linear) using the predefined calibration curves that are suitable for Hach reagents and standard methods, however, the Hach DR5000 spectrophotometer allows the user to re-calibrate the device for any other reagents/ methods by specifying the measurement parameters; i.e. wavelength, relational absorbance values of a series of known concentration, and thus enabling the device to generate the formula that links absorbance values to concentration (the actual work the device does behind the scene) and thus to calculate the concentration of any given sample within the valid</p>	<p>1. الإجراءات الطارئة لمعالجة الحروق الناجمة عن التماس مع السطوح الحارة</p> <p>أ.</p>

formula range.

This process is possible on DR5000 spectrophotometer by using the **User Programs** which enables the user to:

- Programming of user-created procedures.
The user must define or determine the program sequences, calculation formulas, measurement wavelengths, factors, measuring range limits, etc.
- Specifically modified tests.
- Creation of a specific selection of methods and tests.

There are 3 possible methods to calibrate the device

- 1- Entering values; i.e. filling the Concentration- Absorbance table using data provided by reagent manufacturer.
- 2- **Reading Absorbance values of a series of known concentration of standard solution (which is described in this manual)**
- 3- Entering the Concentration/ Absorbance formula provided by the method composer/ reagent manufacturer.

5. Outline of the Process:

- 1- Identifying the parameter which the device will be calibrated to measure.
- 2- Studying the analyses method (as indicated by the reagent manufacturer).
- 3- Preparation of necessary reagents.
- 4- Preparation of enough quantity of standard solutions; this is governed by;
 - a- the measurement method (preparation of the sample, cell volume... etc.)

b- number of values to be used for derivation
of the calibration formula


Remark: *the bigger the number of referral points available for the device to derive the calibration formula, the higher the representation of the calibration curve to the measurement range, and especially when calibration formula is supposed to be polynomial, or 3rd order curve. Usually the curve is linear (straight), within the valid range of measurement (which is governed by the reagents specifications).*

5- Crating the calibration table by entering the concentration values of the standard solutions and then measuring the solutions to determine the corresponding absorbance values. The absorbance values are plotted against the concentrations of the standard solutions and the calibration curve is displayed as a graph.


Note: *In the case of most common parameters, and time-independent, single-wavelength measurements, the calibration curve is a straight, and 4 points (including (0,0) are enough to derivate the representative curve)*

Important: *prepared standard solutions are better to be of regular increment, and equally distributed within the valid range of measurement as much as possible, in order to minimise the effect of extreme collided and values on the curve deviation.*


An exception of the case is when the user would like to create calibration curve for sub-spans of the original span for quality control of the measurements of a limited range of concentration.

Good referral point distribution, points are distributed evenly within the measurement span.

Bad Point Selection, Points are distinct, and no points inside the span (intermediate referral points)




Points are congregated around a certain part of the total valid measurement span





6- Entering the programming mode, and setting up the method (program)


6. Quality Control:




The quality control is conducted by measuring a set of other standard solutions, and thus eliminating the possibility of receiving misleading readings due to the invalidity of the standard solutions used for calibration.

Step	Operation	Remarks(ملاحظات)	التشغيل
1	From the Main Menu, press USER PROGRAMS and then PROGRAM OPTIONS . The Program Options menu will be displayed		
2	Select NEW to program a new user program. The first time PROGRAM OPTIONS is selected, only the NEW option is available. The other options remain inactive (gray) until the first program has been created.		
3	Enter a program number between 949 and 999. The lowest available number appears automatically. The program number will be displayed in the User Program and Favorite Program menu lists (if the method is stored as a favorite). Press OK to confirm. <i>Note: If the program number is already assigned to another user program, a message appears asking whether the existing program should be replaced. Press OK to overwrite an existing program.</i>		



Step	Operation	Remarks(ملاحظات)	التشغيل
4	Use the alphanumeric keypad to enter a program name. The name can be a maximum of 28 characters long. Press BACK to return to the previous screen or press NEXT to continue with the input of the program data.		
5	Select the type of program (Single Wavelength), and press NEXT to continue.		
6	Select the required unit from the list and press NEXT . A user-specific unit that is not included in this list can be added in the edit program under PROGRAM OPTIONS>EDIT .		

Step	Operation	Remarks(ملاحظات)	التشغيل
7	Select the concentration resolution (number of decimal places) and press NEXT .		
8	Enter the chemical form (the chemical representation of the analysis parameter in the result display) using the alphanumeric keypad. Press NEXT to continue.		
9			
10	A method is calibrated by determining the absorbance values of several standard solutions of known concentration. Select "Read Standards" and press NEXT .		

Step	Operation	Remarks(ملاحظات)	التشغيل
11	<p>1. Press the + symbol and use the alphanumeric keypad to enter the standard concentration. Press OK.</p> <p>2. Press the + symbol again and enter the next standard concentration. Repeat this sequence until all standard concentrations (maximum of 24 solutions) are entered.</p> <p>3. Activate the line with the appropriate concentration and insert the cell with the corresponding standard solution.</p>		
	<p>1. Prepare the zero solution and other concentrations according to the intended method (refer to other SOPs if necessary)</p> <p>2. Place the zero solution in the cell compartment. Close the compartment. Press ZERO.</p> <p>3. Place the first standard solution in the cell compartment and close the cell compartment. Press READ.</p> <p>4. Place the second standard solution in the cell compartment and close the cell</p>	<p><i>The TIMER icon shown in the display helps to ensure, when necessary, that the steps of an analysis are correctly timed (e.g. reaction times, wait times, etc., can be exactly specified). When the specified time has elapsed, an acoustic signal is emitted. The use of the timer has no influence on the measurement program.</i></p>	

Step	Operation	Remarks(ملاحظات)	التشغيل
	<p>compartment. Press READ. Repeat this sequence until all the standard solutions are measured (maximum of 24 solutions).</p> <p>5. The entered and measured data are displayed in the Table. To delete a standard concentration, highlight the appropriate line and press DELETE.</p>		
	<p>6. When the data are entered and the measurements are complete, press GRAPH to display the curve that results from plotting the entered data.</p> <p>7. The linear curve corresponds to the standard setting. Press NEXT CURVE to display the polynomial 2nd order curve. Press NEXT CURVE again to display the polynomial 3rd order curve.</p>	 <p>Examples of Good and bad linear curves:</p>  <p>Bad linear graph, curve fit $r^2 < 0.9800$</p>  <p>Good Graph Curve fit $r^2 > 0.9900$</p>	

Step	Operation	Remarks(ملاحظات)	التشغيل
	<p>Upper and Lower Measuring Range Limits</p> <p>Enter the Upper and Lower Limits of the measuring range. An error message is displayed if a reading is above the upper limit or below the lower limit.</p> <ol style="list-style-type: none"> 1. Activate the appropriate line in the overview of the program data and press EDIT. 2. Select On. Press 0.000 to enter the measuring range limit. Press OK to confirm. 		
	<p>Timer Intervals</p> <ol style="list-style-type: none"> 1. Enter the time intervals for the timers. This function can be used to enter abbreviations for work steps and defined intervals for up to four timers. 2. Select the appropriate program in the program data overview screen and press EDIT. 3. Check the appropriate timer(s). Press the TIMER key to select from a list of names that describe the corresponding step. Press 00:00 to enter the times for each timer. 		

Step	Operation	Remarks(ملاحظات)	التشغيل
	<p>Store the program: Review all parameters, when it is deemed necessary to change any parameter, select it by pointing and press EDIT After everything is set press Store to store the program in DR5000 internal memory</p>		
	<p>Formula checkup To check the reagent blank interference you may revise the formula. 1. press calibration 2. select Enter formula 3. press Next (a) represents the affect of the reagent, if you're satisfied about the value, you may press cancel twice and exit programming mode, otherwise, you need to repeat zero concentration solution</p>		

Standard Operation Procedure (SOP) For Sampling

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_____ DFEA

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1. Sample collection

Purpose of the environmental factory inspection is to monitor the pollution sources such as factories and industrial sectors for management of pollution sources. Sampling is the first step in water quality monitoring. Care must always be taken to ensure obtaining a sample that is truly representative. Further, the quality of the sample must be maintained from the time of collection to the time of analysis. If the sample is not representative or if the sample has changed in chemical composition between sampling and analysis the efforts in its analysis become futile.

2. Necessary Equipment and Supplies:

- 1) Sampling vehicle.
- 2) Documents for recording data.
 - ① Manifest for sample receipt,
 - ② Factory information sheet,
 - ③ Others
- 1) Sampling tools
 - a) Sample container/bottle (select materials and capacity depend on the parameters analyzed) (see “Sampling Guide”),
 - b) Sampler/Dipper,
 - c) Bucket
 - d) Gloves,
 - e) Rope,
 - f) Icebox with refrigerant/cooling agent,
 - g) Writing implements, Marker,
 - h) Others
- 2) On-site measurement equipment (as needed)
 - a) Thermometer,
 - b) pH meter,
 - c) EC meter,
 - d) DO meter,
 - e) Others
- 3) Chemicals and tools for sample preservation
 - a) Chemicals
 - (ア) Nitric acid (HNO₃),
 - (イ) Sulfuric acid (H₂SO₄),
 - (ウ) Hydrochloric acid (HCl),
 - (エ) Sodium hydroxide (NaOH),
 - (オ) Others
 - b) Dropping pipet
 - c) pH test paper
 - d) Others

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Step (الخطوة)	Operation	Remarks (ملاحظات)	التشغيل
1	<p>Check sampling vehicle, sampling tools and field measurement by using checklist.</p> <p>The works mentioned above should be completed one day before sampling.</p>		
2	<p>Calibrate pH meter in the site before measuring if it is possible. Otherwise, calibrate it before leaving the laboratory. And do not disconnect the electrode from the instrument.</p> <p>Check EC meter if it is calibrated since less than a month. Otherwise, recalibrate it.</p> <p>Calibrate D.O meter in the site before measuring.</p>		
2	<p>Collect water samples at safe place.</p> <p>The sample bottle should be rinsed three times with portions of sample before being filled. However, this does not apply for bacteriological samples.</p>		

Step (الخطوة)	Operation	Remarks (ملاحظات)	التشغيل
3	Labeling the samples collected. Contents of label should include following contents. 1) Date and time of sampling 2) Sample code 3) Others (name of river, lake, factory etc.)	Example: 2005-10-03 09:30 RAK-I-001 Sugar Factory (I: industrial wastewater R: river L: lake R: reservoir Sp: spring Se: sewage G: groundwater)	
4	Preserve the samples collected following the "Suggested Preservation Methods" and put the sample into icebox.		
5	Conduct field measurement and observation by using "Field Record"		
6	Deliver the samples to laboratory staff.		

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□ Sampling Guide/Sample Preservation

Parameter	Analytical Method Adopted	Container	Preservation	Minimum Sample Volume (Preferred Volume)	Holding Time	Container No.
Oil and Grease	Solvent extraction/Infrared absorptiometry	G, wide-mouth (Washed with soap, rinsed with water, and finally rinsed with solvent to remove any residues)	H ₂ SO ₄ (1:1) or HCl (1:1) to pH below 2. (Don't overfill the sample container and don't subdivide the sample in the laboratory. Refrigerate)	1,000 mL (1,000 mL)	28 days	1
Orthophosphate (PO ₄ ³⁻)	Spectral photometric	P or G (Container should have been cleaned with 1:1 HCl and rinsed with Deionized water. Don't use commercial detergents containing phosphate for cleaning glassware used in phosphate analysis)	Analyze samples immediately for best results. If prompt analysis is not possible, preserve samples by filtering immediately and add 5 mL of chloroform per 1 L of sample, and storing at 4°C. The sample should be at room temperature before analysis.	100 mL	48 hrs	2
Nitrogen, Ammonia (NH ₃ -N)	Spectral photometric	P or G	Refrigerate at 4°C for samples to be analyzed within 24 hr. H ₂ SO ₄ (1:1) or HCl to pH below 2. (Analyze as soon as possible, refrigerate. If chlorine is present, add one drop of 0.1 N Sodium Thiosulfate for each 0.3 mg/L Cl ₂ in a 1-liter sample. Refrigerate. Warm samples to room temperature and neutralize with 5 N Sodium Hydroxide (NaOH) before analysis.	500 mL	7 days (28 days: stored sample)	3

Parameter	Analytical Method Adopted	Container	Preservation	Minimum Sample Volume (Preferred Volume)	Holding Time	Container No.
			Correct the test result for volume additions.)			
Surfactants	Spectral photometric	P or G	Store at 4°C or less (refrigerate) Warm to room temperature before testing	(500 mL)	24 hrs	2
Chromium, Total (Cr-T)	Spectral photometric	P or G (acid-washed container)	HNO ₃ to pH below 2 (Approximately 2 mL per liter of the acid. Store preserved samples at room temperature up to six months. Adjust the pH to about 4 with 5.0 N Sodium Hydroxide (NaOH) before analysis. Correct the test result for volume additions)	(300 mL)	6 months	4
Chromium, Hexavalent (VI)	Spectral photometric	P or G (Rinsed with 1:1 HNO ₃ , Refrigerate)	Store at 4°C in neutrality up to 24 hours. Must be analyzed within 24 hours.	300 mL	24 hrs	2
Sulfide (S ²⁻)	Ion selective electrode	P or G	Fill completely and cap tightly. Avoid excessive agitation or prolonged exposure to air. Add 4 drops of 2N Zinc Acetate ((Zn(CH ₃ COO) ₂ · 2H ₂ O)) solution in the sample bottle per 100 mL sample before sampling. And then adjust pH to above 9 by adding Sodium Hydroxide. Refrigerate	100 mL	7 days	5

Parameter	Analytical Method Adopted	Container	Preservation	Minimum Sample Volume (Preferred Volume)	Holding Time	Container No.
Nitrate (NO ₃ ⁻)	Ion selective electrode	P or G	H ₂ SO ₄ (1:1) to pH below 2, Refrigerate	200 mL	48 hrs	3
Nitrite (NO ₂ ⁻)	Spectral photometric	O or G	Analyze as soon as possible. Add 1 mL of chloroform per 1 L sample. Store 4°C in dark place.	100 mL	48 hrs	
Chloride (Cl ⁻)	Ion selective electrode	P or G	No special preservation	100 mL	28 days	2
Fluoride (F ⁻)	Ion selective electrode	P Glass bottles are satisfactory if previously they have not contained high-fluoride solutions. Always rinse bottle with a portion of sample.	No special preservation	300 mL	28 days	2
Cyanide (CN ⁻)	Ion selective electrode	P (Amber) or G (Amber)	NaOH to pH above 12, Refrigerate	500 mL		6
pH	Ion selective electrode	(On site measurement)				
Electric conductivity (EC)	Ion selective electrode	(On site measurement)				
Suspended solid	Filtrate weight					
Settleable solid	Still standing	P or G Be provided that the material in suspension does not adhere to container walls.	Begin analysis as soon as possible because of the impracticality of preserving the sample. Refrigerate sample at 4°C up to the time of analysis to minimize microbiological decomposition of solids. Transportation and short-term storage of sample will not normally affect the results of the test.		Preferably do not hold samples more than 24 hr. In no case hold sample more than 7 days.	

Parameter	Analytical Method Adopted	Container	Preservation	Minimum Sample Volume (Preferred Volume)	Holding Time	Container No.
			Bring samples to room temperature before analysis.			
COD _{Cr}	Open reflux method	P or G (Preferably collect samples in glass bottles.)	H ₂ SO ₄ (1:1) to pH below 2 Test unstable samples without delay. If delay before analysis is unavoidable, preserve sample by acidification to pH≤2 using conc. H ₂ SO ₄ . Preferably acidify any sample that cannot be analyzed the same day it is collected. Blend samples containing settleable solids with a homogenizer to permit representative sampling. Make preliminary dilutions for wastes containing a high COD to reduce the error inherent in measuring small sample volume. Refrigerate	100 mL	7 days	3
Total coliform	Membrane filter technique	P or G (Should be properly sterilized, either with alcohol, or using an autoclave)	HCl (1:1) to pH below 2	100 mL	6 hrs	7
Heavy Metals	Atomic absorption photometry	G or P (Preferably collect samples in plastic bottles.)	HNO ₃ to pH below 2	100 mL		8

P = Plastic (polyethylene or equivalent); G = Glass

☐ Sampling plan should be made based on kinds of the parameters analyzed.