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
(3)

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RESOLVING TEST FOR MIRROR CAMERA

- * Canon X-ray Mirror Camera CXM5- 70
- * Focal Distance 90 cm
- * Grid ratio 8:1, 34 lines/cm
- * Flourescent Screen Kyokko PO
- * X-ray Film Kodak PF, 70 mm
- * Exposure Conditions 60 KV,
 - Center 2.5 mAs
 - Circumference 3.2 mAs
- * Circumference: 165 mm point on the diagonal line from the center

Resolving Power

Center	2.0 LP/mm
Circumference	1.6 LP/mm

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PHOTOTIMER AND MANUAL COMPARISON
(Indirect Radiography)

No.	Acrl (cm)	KV	mAs	MS	Photo	Film Density
1	5	100	1.4	8.8	+	
2	5	100	1.6	8.0	-	↑
3	8	100	2.7	17.0	+	
4	8	101	2.5	16.0	-	
5	10	100	4.1	25.0	+	
6	10	100	4.0	20.0	-	
7	12	100	6.3	38.0	+	
8	12	100	6.3	32.0	-	
9	5	120	0.8	6.6	+	↑
10	5	120	0.8	5.0	-	
11	8	120	1.5	12.0	+	
12	8	120	1.6	10.0	-	
13	10	120	2.3	18.0	+	
14	10	120	2.5	16.0	-	↑
15	12	120	3.4	27.0	+	
16	12	120	3.2	20.0	-	
17	5	140	3.6	28.0	+	not good
18	5	130	0.7	5.3	+	
19	5	130	0.8	5.0	-	↑
20	8	130	1.2	9.8	+	
21	8	130	1.2	8.0	-	
22	10	130	1.8	14.0	+	
23	10	130	1.6	10.0	-	↓
24	12	130	2.7	21.0	+	
25	12	130	2.5	16.0	-	↓

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PHOTOTIMER AND MANUAL COMPARISON (2)
(Indirect Radiography)

Patients	KV	mAs	ms	Photo	Thickness
A	120	1.7	10	+	Standard
	122	1.6	10	-	
B	120	2.1	13	+	Standard
	121	2.0	12	-	
C	120	2.2	13	+	Standard
	124	2.0	12	-	
D	120	1.5	9.6	+	Standard
	119	1.6	10	-	
E	120	1.5	9.2	+	Standard
	119	1.6	10	-	
F	130	2.4	19	+	Fat
	129	2.5	16	-	
G	120	1.7	10	+	Standard
	121	1.6	10	-	

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RADIOGRAPHIC CONDITIONS

I. INDIRECT RADIOGRAPHY

- Focal Distance 90 cm
- Grid ratio 8:1, 34 lines/cm
- Flourescent Screen Kyokko PO
- Resolving Power Center 2.0 LP/mm
 Circumference 1.6 LP/mm

1. PHOTOTIMER RADIOGRAPHY

KV	Back-up mAs	Photo density	
110	5.0	0	thin
120	5.0	0	standard
130	5.0	0	fat

2. MANUAL RADIOGRAPHY

KV	mAs	
110	1.6 - 2.0	thin
120	1.6 - 2.5	standard
130	2.5 - 3.2	fat

- * mAs in case of manual radiography = display mAs in case of phototimer radiography.
- * mAs setting range > 1.0 mAs
- * 0.1 mAs ↑, 2 KV ↓; 0.1 mAs ↓, 2 KV ↑

2. DIRECT RADIOGRAPHY

- Focal distance 180 cm
- Grid ratio 8:1, 60 lines/cm
- Intensifying Screen SRO - 180
- Film TMG

1. PHOTOTIMER RADIOGRAPHY

KV	Back-up mAs	Photo density	
120	8.0	-1	thin
130	8.0	-1	standard
130	8.0	-1	fat

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An Invitation

*Department of Health
cordially invites you to the
Inauguration Rites of the
Reference Laboratory of the
Cebu Chest Center
on August 15, 1994
at 9:00 in the morning
at the Don Vicente Solto Memorial
Medical Center Compound
Cebu City*

*Inauguration Ceremony
Reference Laboratory of Cebu Chest Center
August 15, 1994
9:00 - 12:00 A.M.
Vicente Sotto Memorial Medical Center Compound
Cebu City*

PROGRAMME

Part I		
9:00 - 9:45 A.M.	-----	Mass Rev. Fa. Nestor Manlangit Officiating Priest
Part II		
National Anthem		
Philippines	-----	to be led by Mrs. Colita Auza
Japan		
Welcome Address	-----	Dr. Milagros B. Bacus OIC, DIRFO 7
	-----	Gov. Vicente Dela Serna Prov. of Cebu
Introduction of Guests/ Participants	-----	Dr. Enrique Sancho Chief, Cebu Chest Center
Message	-----	Secretary Juan Flavie Dept. of Health
	-----	Mr. Akihiko Hashimoto Resident Representative JICA Philippine Office
	-----	Dr. Masakazu Aoki Director, Research Institute of Tuberculosis - Japan Anti-Tuberculosis Association
Part III		
Cutting of Ribbon & Blessing	-----	Secretary Juan Flavie Gov. Vicente Dela Serna Mr. Akihiko Hashimoto Dr. Manuel G. Roxas

Sponsors:

Dr. Ma. Corazon Teoxon	Mrs. Leticia Canoy
Dr. Milagros Bacus	Dr. Rosario Banahay
Dr. Eusebio Alquizalas	Dr. Ramon Macceren
Dr. Jesus Fernandez	Mr. Noel Sin
Dr. Tomas Fernandez	Dr. Elaine Teleron
Mr. Eiji Iwasaki	Mrs. Colita Auza
Dr. Enrique Sancho	Dr. Cristina Giango

Master of Ceremonies:

Dr. Elaine Teleron & Dr. Ramon Macceren

(Lunch will be serve at the SIH Alumni Cultural Center after the Ceremony)


From Cebu to the World

The Freeman

VOL. XXVIII NO. 29 • TUESDAY, AUGUST 16, 1994 • CEBU CITY, PHILIPPINES • PAY ONLY P3.75 in Metro Cebu • 48 Pages

Philippine Office together with Health Undersecretary Dr. Manuel Roxas (second from right) and provincial board member Emmanuel Pepito (extreme right) cut the ceremonial ribbon during the inauguration of the Reference Laboratory of Cebu Chest Center. (FREEMAN photo by Rudolph Alcazaren)

Sun-Star

DAILY

NO. 1 IN READERSHIP. IN CEBU AND, AMONG COMMUNITY PAPERS, THROUGHOUT THE COUNTRY.

CEBU CITY, PHILIPPINES • 48 PAGES • P4.00

TB REFERENCE LAB.

Akihiko Hashimoto (above, second from left), Jica resident representative to the Philippines, and Health Undersecretary Manuel G. Roxas cut the ribbon, marking the ceremonial opening of the P5 million TB reference laboratory of the Cebu Chest Center. They are assisted by Cebu Provincial Board Member Emmanuel Pepito and Health Regional Director Milagros Bacus. The laboratory (left) is said to be the first in the country. It was built through a Jica grant and will be run for five years by the Japan Anti-TB Association; after which the DOH will take over the facility.

SUN-STAR PHOTO: CYRIL CAMPOREDONDO

DOH -JICA

THE PUBLIC HEALTH DEVELOPMENT PROJECT
DOH - FO 7, OSMEÑA BLVD., CEBU CITY 6000, PHILS.

TELEPHONE NOS. 224-029
224-031

FAX NO. (32) 224-027

Message from Dr. Aoki, Chairman of the Domestic Committee on the Philippines Public Health Development Project, JICA

It is a great pleasure for me to extend a few words of congratulation on the occasion of the inauguration of the Reference Laboratory of the Cebu Chest Center, on behalf of the Domestic Committee of the Project, JICA. Needless to say, the bacteriology is the most basic and important technology in the fight against tuberculosis, as clearly defined in the new national tuberculosis programme of the Philippines.

The quality of the bacteriological examination is a crucial issue so that it could affect the entire programme; If the poor

DOH -JICA

THE PUBLIC HEALTH DEVELOPMENT PROJECT
DOH - FO7, OSMEÑA BLVD , CEBU CITY 6000, PHILS.

TELEPHONE NOS. 224-029
224-031

FAX NO. (32) 224-027

examination produces many false cases to be treated as TB, then it may add to the work load of the staff and may waste precious resources; At the same time, low quality service may lose the true cases, or sources of infection, which would be a threat to the patients themselves and to the community, and more seriously, the people's confidence on the public service may be lost. The newly constructed Reference Laboratory of the Cebu Chest Center is expected to ensure the quality service in the national tuberculosis programme to avoid such a disaster as mentioned above, through the training of the personnel, supervisory support of the local microscopy centres and reference services.

DOH -JICA

THE PUBLIC HEALTH DEVELOPMENT PROJECT
DOH - FO 7, OSMEÑA BLVD., CEBU CITY 6000, PHILS

TELEPHONE NOS. 224-029
224-031

FAX NO. (32) 224-027

Now that the new and well equipped facility is available, it is our duty to make best use of it. I sincerely hope that the TB control programme should make a great progress through the use of this laboratory under our further cooperation.

SCHEDULE OF TRAINING ON LABORATORY WORKS

August 29 – September 2, 1994 (1st batch)

1. Ms. Leonides Manatad – IPHO Cebu
2. Ms. Yolanda Roset – IPHO Cebu
3. Ms. Alimay Baron – Aloguinsan
4. Mr. Mario Binoya – Bogu II
5. Ms. Mila Yrauda – Oslob
6. Mr. Alexander Lelis – Lapu-lapu City
7. Ms. Nelle Perez – Mandaue City

September 12 – 16, 1994 (2nd batch)

1. Ms. Marilyn Sua-an – IPHO Cebu
2. Ms. Leonida Clavano – Alcantara & Ronda
3. Mr. Emeritus Reyes – Argao
4. Mr. Raul Galano – Medellin
5. Ms. Felicitas Lumapas – Sogod
6. Ms. Ma. Cecile Evardo – Lapu-lapu City
7. Ms. Jojie Ann Fanlo – Mandaue City
8. Ms. Arleen del Fuero – Mandaue City

September 26 – 30, 1994 (3rd batch)

1. Mr. Felito Villa – IPHO Cebu
2. Ms. Ester Casas – Bogu I
3. Ms. Lourdes Cabahug – Carmen
4. Ms. Lilia C. Tura – Catmon & Tabogon
5. Ms. Eufemia Villahermosa – Dalaguete I & II
6. Ms. Anabelle Buot – Lilo-an
7. Mr. Julfi Tanduyan – Lapu-lapu City
8. Ms. Maria Edna Cortes – Mandaue City

October 3 – 7, 1994 (4th batch)

1. Ms. Yvonne Cariñosa – IPHO Cebu
2. Ms. Rhoda Cañete – San Remegio
3. Ms. Lilia Gonzales – Badian
4. Mr. Aaron Alvarez – Lapu-lapu City
5. Mr. Rogelio Dutallas – Lapu-lapu City
6. Ms. Teresita Ortalla – Mandaue City
7. Mr. Kerwin Naqua – Mandaue City
8. Mr. Ephraim Gonzaga – Barili District Hospital

NTP Refresher Course on Laboratory Works
Reference Laboratory of Cebu Chest Center
August 29 - September 2, 1994

P R O G R A M M E

<i>Invocation</i>	-----	<i>Ms. Nelle Perez</i>
<i>Philippine National Anthem</i>	-----	<i>Ms. Leonides Manatad</i>
<i>Opening Remarks</i>	-----	<i>Dr. Ramon Maceren</i>
<i>Message</i>	-----	<i>Dr. Milagros B. Bacus</i> <i>OIC, DIRFO 7</i>
<i>Introduction of Participants</i>	-----	<i>Ms. Colita Auza</i>
<i>Mechanics of the Seminar</i>	-----	<i>Dr. Elaine Teleron</i>
<i>Moderator</i>	-----	<i>Ms. Colita Auza</i>

SCHEDULE OF ACTIVITIES

Day	AM 8:00 – 12:00 (15' Tea Break)	PM 1:30 – 5:00 (15' Tea Break)
1	8:00 – 8:30 Registration 8:30 – 8:45 Opening Ceremony Course Orientation 8:45 – 9:30 Quiz (Pre-Test) 9:45 –12:00 (Pr) What are the difficulties encountered in your laboratory?	1:30 – 5:00 (P) Laboratory Preparation – Laboratory set up – Reagents preparation – Microscope handling – Evaluatin of Ziehl-Neelsen staining techniques (P) Ziehl-Neelsen staining
2	8:00 – 8:55 (L) General concept on TB and TB control 9:00 – 9:55 (L) Aim of TB Lab. services 10:10 –12:00 (L) Specimen Collection (L) Safety precautions & disposal	1:00 – 5:00 (L) Ziehl-Neelsen method (P) Chronological studies (1) (P) Sputum Collection
3	8:00 –12:00 (L) NTP in the Philippines (L) Recording/reporting system (L) Logistics	1:30 – 5:00 (P) Chronological studies (2) (P) Smear slide reading
4	8:00 –12:00 (P) Chronological studies (3) (P) Distribution of AFB in the sputum	1:30 – 5.00 (L) Quality check for direct smear examination (P) Quality check for direct smear examination – Smear slide checking – Microscopy reading check
5	8:00 –10:45 (P) Chronological studies (4) 11:30 –12.00 Quiz (Post-test)	1:30 – 2:30 Course Evaluation and Summary Closing Ceremony

(L) Lecture, (P) Practice, (Pr) Presentation

LIST OF REAGENTS AND EQUIPMENTS NEEDED FOR MICROSCOPY TRAINING

Reagents	Minimum quantity
1. Fuchsin basic	1.5 g/group
2. Methylene blue	0.3 g/group
3. Sulfuric acid (H ₂ SO ₄)	100 ml/group
4. Phenol	15 ml/group
5. Ethanol	130 ml/group
6. Xylene	400 ml/group
7. Immersion oil	5 ml/person
8. Distilled water	1000 ml/group

LIST OF REAGENTS & EQUIPMENTS NEEDED FOR MICROSCOPY TRAINING

Glasswares	Minimum quantity
1. Flask	
(500 ml)	3/group
(1000 ml)	1/group
2. Cylinder	
(50 ml)	1/group
(200 ml)	2/group
3. Pipette (10 ml)	1/group
4. Reagent bottle	
(Brown, 300 ml)	2/group
(Transparent, 500 ml)	1/group
5. Funnel	1/group
6. Glass slide	20/person
7. Plate (Petri dish)	1/person
8. Staining jar	1/person

LIST OF REAGENTS & EQUIPMENTS NEEDED FOR MICROSCOPY TRAINING

Equipments	Minimum quantity
1. Filter paper	5 sheets/person
2. Forceps	1/person
3. Coco midribs	20/person
4. Loop holder	1/person
5. Nicromwire	1/group
6. Tally counter	1/person
7. Washing bottle	2/clean bench
8. Spirit cotton	1/room
9. 5% phenol (disinfectant sol.)	1/room
10. Slide rack	1/person
11. Staining bridge	2/group
12. Water bath with tripod	1/room
13. Pliers	2/group
14. Microscope	1/person
15. Gauze/lenspaper	few/person
16. Lens cleaner set (Xylene/petrorium benzine, Ethanol)	1/group
17. Smear slide tray (mappe)	2/person
18. Safety pipetter	1/room
19. Staining solution dropper set	1/group
20. Vinyl tape	1/group
21. Ruler	2/group

National Tuberculosis Control Program
REFRESHER TRAINING ON LABORATORY WORKS

I. PURPOSE:

This training is carried out to address the needs of Medical Technologists for a refresher training on Laboratory works. The National Tuberculosis Control Program of the Department of Health adheres to provide good services to the populace through early detection of tuberculosis cases in the community. As a means to achieve this goal, it uses a simple, affordable and accurate diagnostic tool, which is the Direct Smear Examination. With this, it is imperative that Medical Technologists that manned the microscopy centers and the key personnel in TB diagnosis should be adept with knowledge and skills in TB Control and laboratory techniques with emphasis on sputum microscopy.

This refresher program accomplished in a series of batches enables participants from the DOH-JICA Intensive Service Areas (ISA) and Cebu, IPHO validators to implement a unified and strengthened NTP in Cebu Province.

II. OBJECTIVES:

A. General Objectives:

- 1) To provide inputs on TB Control and laboratory techniques on Direct Smear Examination
- 2) To upgrade their knowledge and skills in sputum microscopy through practice examination.

B. Specific Objectives:

By the end of the training period, the participants are expected to be able to:

- 1) Know the importance of direct smear examination in the TB control programme
- 2) Gain insights on the techniques of sputum collection, smearing, fixation, staining and microscopy reading
- 3) Manage and evaluate smear microscopy laboratory
- 4) Adopt a system of recording/reporting, logistics and referral between laboratories/radiographic facilities

III. METHODOLOGY:

Laboratory Practice
Lecture Discussion
Presentation
Pre and post test
Demonstraion/Return Demonstration

IV. OPERATING DETAILS:

Duration : August 29 – September 2, 1994 (1st Batch)
September 12 – September 16, 1994 (2nd Batch)
September 26 – September 26, 1994 (3rd Batch)
October 3 – October 7, 1994 (4th Batch)

Venue : Reference Laboratory of Cebu Chest Center

Time : 8:00 a.m. – 5:00 p.m.

Participants : See attached list of Medical Technologists

Facilitators:

August 29 – September 2, 1994

Dept. of Health & Integrated Regional Field Office No. 7

Dr. Nora Cruz – Training Officer, TB Control Service, Manila
Ms. Paz Rostrata – Medical Technologist, TBCS, Manila
Dr. Elaine Teleron – Regional TB Medical Coordinator
Ms. Colita Auza – Regional TB Nurse Coordinator
Mr. Benny Loberiza – Medical Technician, Cebu Chest Center
Dr. Enrique Sancho – Chief, Cebu Chest Center
Ms. Leonides Manatad – Medical Technologist, IPHO Cebu

Japan Internatinal Cooperation Agency

Dr. Masashi Suchi – DOH–JICA Chief Adviser
Ms. Akiko Fujiki – Short–term Expert on Laboratory Works

Funding : DOH–JICA Project

V. REQUIREMENTS:

Complete Attendance
Active participation

CURRICULUM AND OBJECTIVES

LECTURE :

1) TB bacteriology and TB laboratory works

To review basic TB bacteriology and TB laboratory works

- Ziehl-Neelsen method (Mr. Loberiza)
- Specimen collection (Mr. Loberiza)
- Safety precautions and disposal (Ms. Rostrata, Mr. Loberiza)
- Quality check for direct smear examination (Ms. Fujiki, Mr. Loberiza)

2) TB Control

To review general concept of TB control and role of bacteriology examination

- General concept on TB and TB Control (Dr. Teleron, Dr. Sancho)
- Aim of TB Laboratory Services (Dr. Cruz, Ms. Rostrata, Mr. Loberiza)
- NTP in the Philippines (Dr. Cruz, Dr. Teleron, Dr. Sancho)
- Recording/reporting system (Dr. Cruz, Ms. Auza, Dr. Sancho)
- Logistics (Dr. Cruz, Ms. Rostrata, Mr. Loberiza, Ms. Auza, Dr. Sancho)

PRACTICE :

1) Laboratory preparation (Mr. Loberiza, Ms. Manatad)

- a) To learn how to fix and prepare equipments and reagents needed for microscopy laboratory
- b) To gain the knowledge and learn how to manipulate and maintain a microscope
 - Laboratory set-up
 - Reagent preparation
 - Microscope handling
 - Evaluation of Ziehl-Neelsen staining techniques

2) Ziehl-Neelsen staining and Direct smear examination (Mr. Loberiza, Ms. Manatad)

To learn the techniques of detection of AFB by Ziehl-Neelsen method and to observe the AFB in the sputum specimens

3) Smear slide reading (Mr. Loberiza, Ms. Manatad)

To know inter/intra variations among the persons who read the results of direct smear

4) Chronological studies (Mr. Loberiza, Ms. Manatad)

To observe the chronological changes of the sputum specimens and the AFB after storage at room temperature

5) Distribution of AFB in the sputum (Mr. Loberiza, Ms. Manatad)

To observe the variety of number of AFB in the parts of sputum specimen by direct smear examination.

- 6) Quality check for direct smear examination (Mr. Loberiza, Ms. Manatad)
To learn how to check and evaluate stained smear slides

PRESENTATION :

What are the difficulties encountered in your laboratory? (All Facilitators)

To know the actual encountered problems in different laboratories and to share the knowledge and experiences for the solution of the problems

THE GENERAL RULES OF LABORATORY CONDUCT

1. Make sure that all gas, water, and electrical appliances at your table are turned off before you leave.
2. Eating, drinking or smoking is prohibited in the laboratory at all times.
3. If sputum specimen or a living culture of microorganisms is spilled, notify the instructor immediately.
4. If you are injured in any way, notify the instructor immediately.
5. Wash your hands at the conclusion of each laboratory period.
6. Do support the microscope carefully, when moving it.
7. Be sure you wipe off the immersion oil before putting your microscope away.
8. If you destroyed the laboratory equipment or material, notify the instructor immediately.

You are strictly requested to keep above mentioned rules in the practice room.

REAGENTS PREPARATION FOR ZIEHL-NEELSEN STAIN

(1) Carbol Fuchsin (330 ml/group)

1) Stock Solution

- * Basic Fuchsin 1.5 g
- * Ethyl alcohol 30 ml

2) 5% Phenol Solution

- * Phenol, melted crystals 15 ml
- * Distilled water 285 ml

3) Carbol Fuchsin Solution

- * Stock solution 30 ml
- * 5% Phenol solution 300 ml

(2) 20% Sulfuric acid (H_2SO_4) solution (500 ml/group)

- * Sulfuric acid 100 ml
- * Distilled water 400 ml

(3) 0.1% Methylene Blue (300 ml/group)

- * Methylene blue 0.3 g
- * Distilled water 300 ml

(Memo)

ZIEHL-NEELSEN STAIN

Specimen : Fresh sputum containing AFB

Reagent : Carbol fuchsin solution
 20% H₂SO₄
 0.1% Methylene blue solution

Procedure : (1) Smear Preparation
 1) Smearing using a coco midrib and wire loop, 10 each
 2) Drying
 3) Fixation
 (2) Staining by Ziehl-Neelsen method
 1) Carbol fuchsin solution
 2) Rinse off
 3) 20% H₂SO₄
 4) Rinse off
 5) Methylene blue
 6) Rinse off and blott to dry
 (3) Reading
 1) Select 2 smear slides, one from each group
 2) Scan 2 horizontal lines or 3 vertical lines
 3) Count and record the number of AFB in each of 10 fields
 4) Sketch the AFB

(Memo)

VF	No. of AFB		VF	No. of AFB		Sketch the AFB and background of smear
	loop	coco		loop	coco	
1			6			
2			7			
3			8			
4			9			
5			10			

Total No. of AFB in 2 horizontal lines/3 vertical lines by wire loop _____
 by coco midrib _____

CHRONOLOGICAL STUDIES ON STORED SPUTUM

Specimen : Fresh sputum containing AFB

- Procedure :
- 1) Spread sputum specimen on a slide with a coco midrib stick
 - 2) Fix and stain by Ziehl-Neelsen method
 - 3) Seal the sputum specimen and store at room temperature
 - 4) Observe macroscopically and microscopically the following points:
 - Appearance of sputum specimen
 - Number of tubercle bacilli
 - Stained condition and forms of AFB
 - Background of the stained film
 - Anything that you observe and recognize
 - 5) Repeat the above procedures on the following days:
 - 1st day after 1st preparation
 - 2nd day after 1st preparation
 - 3rd day after 1st preparation

(Memo)

Duration	Appearance of sputum colour viscosity	No. of AFB		Stained and form condition of AFB	Background
		loop	coco		
Day of preparation					
1st day					
2nd day					
3rd day					

DISTRIBUTION OF TUBERCLE BACILLI IN SPUTUM

Specimen : Fresh sputum containing AFB

Procedure : 1) Pour the sputum into the plate
2) Observe the properties of the sputum specimen and sketch them
3) Fish the 5 different specks of sputum (with a coco midrib stick)
4) Make the direct smears from each speck and stain them by Ziehl-Neelsen method
5) Count the number of AFB in each stained smear
6) Compare the properties of sputum in each area and the number of AFB seen in each smear

(Memo)

Specimen	Macroscopic Appearance	No. of AFB	Macroscopic Whole Appearance
1			
2			
3			
4			
5			

SMEAR SLIDE READING

Specimen : 20 stained smear slides

- Procedure :
- 1) Read the 20 smear slides
 - 2) Record the results by following scale
 - (-) No AFB seen in whole smear
 - (+) 1 - 5 AFB seen in whole smear
 - (++) 6 - 24 AFB seen in each field
 - (+++) 25 or more AFB seen in each field
 - 3) Analyse the inter and intra variation of smear reading in statistical way

(Memo)

No.	Your Result	No.	Your Result
1		11	
2		12	
3		13	
4		14	
5		15	
6		16	
7		17	
8		18	
9		19	
10		20	

Result of standard reader	Your Result					Total
	-	+	++	+++	++++	
-						
+						
++						
+++						
++++						
Total						

Name : _____

Encircle the correct answer to each question.

1. Tubercle bacillus was discovered in 1882 by
 - a. R. Koch
 - b. L. Pasteur
 - c. P. Erlich
 - d. S. A. Waksman
 - e. L.A. Calmette
2. What is the color of AFB when stained by Ziehl–Neelsen method?
 - a. red
 - b. blue
 - c. yellow
 - d. black
 - e. green
3. The most adequate magnification for ordinary microscoping of tubercle bacilli in sputum smear is
 - a. 40 x 10
 - b. 100 x 10
 - c. 100 x 15
 - d. 10 x 10
 - e. 20 x 5
4. The term "acid–fast" as applied to mycobacteria means that :
 - a. The organisms produce acid from sugar very rapidly
 - b. Once stained, the organisms resist decolorization with acids
 - c. The organisms are resistant for disinfection by acid.
 - d. None of the above
 - e. All of the above
5. How many visual fields should be examined by x1000 magnification before reporting an acid–fast smear negative?
 - a. 300 VF
 - b. 100 VF
 - c. 40 – 50 VF
 - d. 1,000 VF
 - e. 150 VF
6. Which reagent is not used for Ziehl–Neelsen stain?
 - a. Methylene blue
 - b. Auramine
 - c. Fuchsin
 - d. Phenol
 - e. Hydrochloric acid
7. Which manner is proper during sputum smear examination?
 - a. Reporting red stained rods as "tubercle bacilli" after acid–fast microscopy
 - b. Heating the slide flood with Ziehl's solution to boiling
 - c. Staining sputum smears from different patients together in a stained jar.
 - d. Flaming a wire loop after spreading sputum on a slide for sterilization
 - e. None of the above
8. Which route is most common in tuberculosis infection?
 - a. swallowing in digestive tract
 - b. Inhalation in respiratory tract
 - c. Invasion from wound in skin
 - d. Sticking in skin
 - e. Splashing into eyes

9. The "sputum" that is most productive of positive smear results is described as:
- a. mucous
 - b. clear and watery
 - c. mucopurulent or purulent
 - d. saliva
 - e. hemoptysis

10. Which one is false?

- a. Ordinary surgical mask will protect against tuberculosis infection.
- b. Using a biological safety cabinet is one of the important part of laboratory safety programme.
- c. Microscopy provides a rapid method to establish a precise diagnosis of tuberculosis
- d. Acid - fast bacilli in smears lose their acid-fastness by U.V. light and excess heat
- e. Carbol fuchsin and fluorochrome staining methods are mostly used to demonstrate the presence of mycobacteria in smears.

1. How long have you engaged in the practical work of laboratory?

_____ years

2. How many specimens do you receive for TB smear examination?

_____ specimens/day/week/month

3. Draw the actual size of smear and indicate how you scan on the smear.

4. What is your reporting scale?

-- _____ ++ _____

± _____ +++ _____

+ _____ +++++ _____

EVALUATION OF PRE-POST TEST FOR FIRST BATCH, 1994

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B: Before training : Average 65.7

A: After training : Average 91.4

EVALUATION FOR REFRESHER TRAINING ON LABORATORY WORKS
(August 29 – October 7, 1994)

1.) Course

Duration	:	Just right	Too short	Too long
Conduction	:	Good	Fair	Poor
Facilities	:	Good	Fair	Poor
Tea break	:	Good	Fair	Poor
Lunch	:	Good	Fair	Poor
Accommodation	:	Good	Fair	Poor
Allowance	:	Just right	Fair	
Benefit obtained	:	Yes	No	Don't know

2.) Session/Lectures

Time allocation	:	Just right	Too short	Too long
Level	:	Just right	Too high	Too primary
Content	:	Just right	Too broad	Too specific

3.) Reasons to answer "Fair/Poor" or "No/I Don't Know" are :

4.) Topics which you liked most are :

5.) Topics which you want to learn more are :

6.) Topics which are difficult to understand are:

7.) Other comments

EVALUATION FOR REFRESHER TRAINING ON LABORATORY WORKS
(August 29 - September 2)

1) Course

Duration:	Just right	6	Too short	1	Too long	0
Conduction:	Good	7				
Facilities:	Good	7				
Tea break:	Good	4	Fair	3		
Lunch:	Good	7				
Accommodation:	Good	7				
Benefit obtained:	Yes	7				

2) Session/Lectures

Time allocation:	Just right	7				
Level:	Just right	7				
Content:	Just right	5	Too broad	1	Too specific	1

3) Reasons to answer "Fair/poor" or "No/don't know" are:

- Tea break - serve very late in the afternoon

4) Topics which you liked most are:

- The preparation of AFB stain
- The chronological study of sputum specimen
- General concept on TB and TB control
- Aim of TB laboratory services
- Specimen collection
- Safety precautions and disposal
- Logistics
- Direct smear examination
- All the topics
- Differentiating *M. bovis* and *M. tuberculosis* and the staining, semaring of AFB
- Discussion like the one, in which the field or RIU med. tech complain that there reagents are sub-standard and we have found out that it is only in the timing of heating and staining the specific slides.
- The collection of the sputum, the problems meet in the field in a sputum microscopy, the preparation of the reagent, direct smear examination
- Mr. Loberiza is good in expressing the produce.

5) Topics which you want to learn more are:

- None
- General concept on TB and TB control
- Mycobacterium other than TB, the differentiation of each.
- Discuss more on the size and shape of the bacilli after delayed smearing of sputum
- About validation, I want to learn more because I want that I could apply it in the field sufficiently and diligently
- Topic on culture and sensitivity so that we can do in our patient which are resistant to the latest TB drugs

6) Topics which are different to understand are:

- None
- What reporting scale to follow the DOH-JICA on the National standard of the Philippines
- So far, so good because we have a very nice lectures, Mdm. Akiko Fujiki, Dr. Nora Cruz and she explain to asked if ever we could not understand.

7) Other comments

- The foods served during tea break is not so good and sometimes water is not available
- No comment. Thank you.
- Very small time given in scanning smears
- None
- I have no comments about the training, because it is our long dream to have a proper training like this one, and I think this kind of training are very worth to all of us who are working for so long. We thank to our facilitators, Dr. Nora Cruz, Dr. E. Teleron, Paz. Benny and to Dr. Suchi, Mdm. Akiko Fujiki, Mr. Terasaki and to JICA.
- No comments

Closing Programme

Refresher Course on Laboratory Works
September 2, 1994

Invocation	Ms Almay Baron
Pambansang Awit	Ms Lenny Manatad
Impressions From Participants	Ms Nelly Perez
From Facilitators	Mr Benny Loberiza
Message	Dr Ma Corazon Teoxon OIC, TB Control Service
Closing Remarks	Dr. Masashi Suchi Chief Adviser DOH-JICA Project Office

Emcee Ms Yolanda Garces